

MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MULTICAM](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 100050369
- Mfr. No.: HW-TKN-MTKN-ESS-MC
- Size: 1 Person
- Delivery weight: 0.322kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650022

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Micro Trauma Kit NOW! Essentials Supplies Safety Instructions](#)
- [About Us](#)

Micro Trauma Kit NOW! Essentials Supplies

Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Essentials Supplies by Blue Force Gear. This kit is designed for everyday carry and contains essential lifesaving medical supplies. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure the kit is kept in a secure location where it is easily accessible in case of emergencies.
- Regularly check the contents of the kit to ensure all items are present and in good condition.
- Use the kit only for its intended purpose—providing first aid in emergency situations.
- Keep the kit out of reach of children and vulnerable individuals unless under adult supervision.
- Report any unsafe products or accidents to the appropriate authorities.
- Stay informed about product recalls by checking the EU Safety Gate platform.

Specific Safety Precautions for Use

- Always wear trauma gloves when providing first aid to prevent contamination.
- Use the Quickclot Bleeding Control Dressing only as directed on the packaging.
- Ensure that the Mini Compression Dressing is applied correctly to avoid further injury.
- The Hypothermia Blanket should be used to maintain body temperature; ensure it is placed correctly around the individual.
- Use Flat Fold Tape to secure dressings, ensuring it does not restrict blood flow.
- If you are unsure how to use any item in the kit, seek guidance from a qualified medical professional.

Instructions for Installation and Usage

1. Familiarize Yourself with the Kit:

- Open the pouch and review the contents.
- Identify each component and understand its purpose.

2. Accessing Supplies:

- The kit features BLIP pull tabs for easy access.
- Pull the tab with one hand or a single finger from either side to deploy supplies.

3. Using the Trauma Gloves:

- Before providing aid, put on the trauma gloves to protect yourself and the patient.

4. Applying the Quickclot Dressing:

- Remove the dressing from the kit.
- Place it directly on the wound and apply pressure.

5. Using the Mini Compression Dressing:

- Place the dressing over the Quickclot dressing if bleeding persists.
- Secure it in place with Flat Fold Tape.

6. Using the Hypothermia Blanket:

- Unfold the blanket and wrap it around the individual to help maintain body temperature.

7. Securing Dressings:

- Use Flat Fold Tape to secure all dressings in place, ensuring they are tight but not overly restrictive.

8. After Use:

- Dispose of used gloves and dressings properly.
- Restock any used items in the kit as soon as possible.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Ensure that any sharp objects, such as used dressings, are disposed of safely to prevent injury.
- Recycle the packaging materials where possible.

Contact Information for Further Support

For any inquiries regarding the Micro Trauma Kit NOW! Essentials Supplies, please reach out to your local distributor or visit the official Blue Force Gear website for more information.

Thank you for prioritizing safety with the Micro Trauma Kit NOW! Essentials Supplies. Stay prepared and safe!

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk