

MICRO TRAUMA KIT NOW! PRO SUPPLIES - BLUE FORCE GEAR

MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT BLACK

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components - the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The PRO Kit Includes: Trauma Gloves (1x Pair) Quickclot Combat Gauze LE (1x) SWAT-T Pressure Device (1x) Hypothermia Blanket - Green/Silver (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT BLACK](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 100050373
- Mfr. No.: HW-TKN-MTKN-PRO-BK
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650053

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Micro Trauma Kit NOW! Pro Supplies Safety Instructions](#)
- [About Us](#)

Micro Trauma Kit NOW! Pro Supplies Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Pro Supplies by Blue Force Gear. This guide provides essential safety instructions to ensure the safe use of this trauma kit. Please read this manual carefully before use to understand how to handle the kit properly and avoid potential hazards.

General Safety Guidelines

- Ensure that the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure that all items are present and in good condition.
- Replace any expired or damaged items immediately.
- Keep the kit out of reach of children unless under adult supervision.
- Familiarize yourself with the contents of the kit and their intended use.
- In case of an emergency, call for professional medical assistance immediately after using the kit.

Specific Safety Precautions for Use

- **Trauma Gloves:** Wear gloves when providing first aid to prevent contamination and protect yourself from bloodborne pathogens.
- **Quickclot Combat Gauze:** Use only as directed for severe bleeding. Apply direct pressure after placement.
- **SWATT Pressure Device:** Ensure proper application to control bleeding. Follow instructions carefully.
- **Hypothermia Blanket:** Use to maintain body heat in case of shock or exposure to cold. Ensure the patient is properly wrapped.
- **Age Restrictions:** This kit is designed for use by adults and trained professionals. Use caution when allowing minors to handle any items.

Instructions for Installation and Usage

1. Accessing the Kit:

- The Micro TKN can be mounted on a belt or MOLLE system. Ensure the mounting is secure before use.
- Use the BLIP pull tabs for quick access to the supplies with one hand.

2. Using the Supplies:

- **Trauma Gloves:**
 - Open the glove package and put on gloves before any first aid procedures.
- **Quickclot Combat Gauze:**
 - Remove from packaging, apply directly to the wound, and apply pressure.
- **SWATT Pressure Device:**
 - Unroll and wrap around the injured area, securing it tightly to control bleeding.
- **Hypothermia Blanket:**
 - Unfold and wrap around the individual needing warmth, ensuring full coverage.

3. After Use:

- Dispose of used items according to local regulations.
- Replenish supplies as needed to maintain kit readiness.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations for hazardous waste.
- Contact your local waste management authority for guidance on proper disposal of medical items.
- Do not throw away gloves or other contaminated items in regular trash.

Contact Information for Further Support

For any questions or concerns regarding the Micro Trauma Kit NOW! Pro Supplies, please refer to the manufacturer's website or contact support through the provided channels.

Conclusion

The Micro Trauma Kit NOW! Pro Supplies is designed to provide essential first aid in emergency situations. By following these safety instructions, you can ensure effective and safe use of the kit for yourself and others. Always prioritize safety and seek professional medical assistance when needed.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk