MICRO TRAUMA KIT NOW! PRO SUPPLIES BELT MOUNT COYOTE

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components – the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The PRO Kit Includes: Trauma Gloves (1x Pair) Quickclot Combat Gauze LE (1x) SWAT-T Pressure Device (1x) Hypothermia Blanket - Green/Silver (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! PRO SUPPLIES BELT MOUNT COYOTE
- Manufacturer: BLUE FORCE GEAR
- Product no.: 100050377
- Mfr. No.: BT-TKN-MTKN-PRO-CB
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650169

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Micro Trauma Kit NOW! Pro Supplies Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Pro Supplies. This compact trauma kit is designed for everyday carry and contains essential lifesaving medical supplies. It is suitable for law enforcement professionals, prepared citizens, and hunters. This guide provides important safety instructions to ensure the safe and effective use of your trauma kit.

General Safety Guidelines

- **Product Safety:** Always ensure that the kit is used according to the instructions provided. The kit is designed to minimize risks associated with medical emergencies.
- **Enhanced Recalls:** Stay informed about any recalls. If you receive a recall notice, follow the instructions to return or replace the product.
- **Online Shopping:** If purchased online, ensure that the seller complies with safety requirements. Always check for product authenticity.
- **Special Consumer Focus:** Keep the kit out of reach of children. This kit contains items that may pose hazards if misused.
- **EU Contact Point:** For any safety inquiries, please refer to the contact information provided with your product.
- **Rapid Alerts:** Regularly check updates on the EU's Safety Gate platform for any alerts related to this product.

Specific Safety Precautions for Use

- **Read Instructions:** Before using the kit, read all instructions carefully.
- **Check Contents:** Regularly check the contents of your kit to ensure all items are present and not expired.
- **Proper Handling:** Handle all medical supplies with care to avoid injury.
- **OneHanded Deployment:** Practice using the BLIP pull tabs to familiarize yourself with onehanded deployment for quick access in emergencies.
- **Environmental Considerations:** Store the kit in a dry, cool place away from direct sunlight to preserve the integrity of the medical supplies.

Instructions for Installation and Usage

1. Belt Mounting:

- Select a suitable location on your belt where the kit can be easily accessed.
- Securely attach the MOLLE or beltmounted pouch using the TenSpeed technology to ensure it does not shift during movement.

2. Accessing Supplies:

- To access the supplies, locate the BLIP pull tabs on either side of the pouch.
- Pull the tab with one hand or a single finger to deploy the insert containing the medical supplies.

3. Using Medical Supplies:

- **Trauma Gloves:** Always wear gloves when administering first aid to prevent contamination.
- **Quickclot Combat Gauze:** Follow the instructions on the packaging for proper usage.
- **SWATT Pressure Device:** Familiarize yourself with the application process for effective use.
- **Hypothermia Blanket:** Use as needed to maintain body temperature in emergency situations.

Disposal Instructions

- Dispose of any expired or used medical supplies in accordance with local regulations.
- Do not dispose of the kit itself in regular household waste. Check for local recycling programs or hazardous waste disposal options.

Contact Information for Further Support

For any questions or concerns regarding the Micro Trauma Kit NOW! Pro Supplies, refer to the contact information provided with your product for assistance.

Thank you for prioritizing safety and preparedness with the Micro Trauma Kit NOW! Pro Supplies. Always be ready to respond to emergencies effectively!

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