

MICRO TRAUMA KIT NOW! PRO SUPPLIES - BLUE FORCE GEAR

MICRO TRAUMA KIT NOW! PRO SUPPLIES BELT MOUNT MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components - the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The PRO Kit Includes: Trauma Gloves (1x Pair) Quickclot Combat Gauze LE (1x) SWAT-T Pressure Device (1x) Hypothermia Blanket - Green/Silver (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! PRO SUPPLIES BELT MOUNT MULTICAM](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 100050378
- Mfr. No.: BT-TKN-MTKN-PRO-MC
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 114mm
- Shipping length: 178mm
- UPC: 810073650176

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Micro Trauma Kit NOW! Pro Supplies Safety Instructions](#)
- [About Us](#)

Micro Trauma Kit NOW! Pro Supplies Safety Instructions

Introduction

Thank you for choosing the Micro Trauma Kit NOW! Pro Supplies from Blue Force Gear. This kit is designed to provide essential lifesaving medical supplies for everyday carry. It is important to follow the safety instructions and guidelines provided in this manual to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and to ensure that all items are in good condition.
- Keep the kit out of reach of children unless under adult supervision.
- Use the kit only for its intended purpose—providing first aid in emergencies.
- In case of any doubt about the use of any item in the kit, consult a medical professional.

Specific Safety Precautions for Use

- Trauma Gloves: Always wear gloves when treating injuries to prevent contamination and protect yourself from bloodborne pathogens.
- Quickclot Combat Gauze: Follow the instructions carefully for application to control bleeding. Do not use if the packaging is damaged.
- SWATT Pressure Device: Ensure proper application according to the included instructions. Do not apply too tightly to avoid cutting off circulation.
- Hypothermia Blanket: Use the blanket to maintain body temperature in cold situations. Ensure that the patient is not allergic to the materials used in the blanket.

Instructions for Installation and Usage

1. Accessing the Kit:

- The Micro TKN is designed for onehanded access. Use the BLIP pull tabs to deploy the contents quickly.
- Practice accessing the kit in a safe environment to familiarize yourself with the process.

2. Using the Trauma Gloves:

- Open the glove pouch and put on the gloves before attending to any injuries.
- Dispose of the gloves properly after use to avoid contamination.

3. Applying Quickclot Combat Gauze:

- Open the package and apply the gauze directly to the wound.
- Apply pressure to the wound for several minutes to help control bleeding.

4. Using the SWATT Pressure Device:

- Unroll the SWATT device and place it over the injury.
- Wrap it around the limb and secure it in place, ensuring it is tight enough to control bleeding but not so tight that it cuts off circulation.

5. Using the Hypothermia Blanket:

- Unwrap the blanket and cover the patient to maintain their body heat.
- Ensure that the shiny side is facing the patient for maximum insulation.

Disposal Instructions

- Dispose of used medical supplies, such as gloves and gauze, in a sealed plastic bag to prevent contamination.
- Follow local regulations for the disposal of medical waste.
- Ensure that the kit itself is disposed of properly if it becomes damaged or is no longer needed.

Contact Information for Further Support

For any inquiries or concerns regarding the Micro Trauma Kit NOW! Pro Supplies, please refer to the official website or customer support resources provided by Blue Force Gear.

By adhering to these safety instructions, you can ensure the effective and safe use of your Micro Trauma Kit NOW! Pro Supplies. Remember, in any emergency, seeking professional medical assistance is always recommended.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk