

# MICRO TRAUMA KIT NOW! ESSENTIALS SUPPLIES - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNT ODG

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components - the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The Essentials Kit Includes: Trauma Gloves (1x Pair) Quickclot Bleeding Control Dressing (1x) Mini Compression Dressing (1x) Hypothermia Blanket - Orange/Silver (1x) Flat Fold Tape (1x)



## Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNT ODG](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430100743
- Mfr. No.: HW-TKN-MTKN-ESS-OD
- Size: 1 Person
- Delivery weight: 0.222kg
- Shipping height: 51mm
- Shipping width: 152mm
- Shipping length: 229mm
- UPC: 810073650039

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [Micro Trauma Kit NOW! Essentials Supplies Safety Instructions](#)
- [About Us](#)

# Micro Trauma Kit NOW! Essentials Supplies

## Safety Instructions

### Introduction

Thank you for choosing the Micro Trauma Kit NOW! Essentials Supplies by Blue Force Gear. This kit is designed to provide essential lifesaving medical supplies for everyday carry. It is important to follow the safety instructions outlined in this guide to ensure safe and effective use of the product.

### General Safety Guidelines

- Ensure that the kit is kept out of reach of children and vulnerable individuals.
- Regularly check the contents of the kit to ensure all supplies are intact and within their expiration dates.
- Use the kit only for its intended purpose: to provide first aid in emergency situations.
- Report any unsafe products or accidents to the appropriate authorities.
- Stay informed about any recalls or safety notices related to this product via the EU's Safety Gate platform.

### Specific Safety Precautions for Use

- Always wear trauma gloves when handling any medical supplies to prevent contamination.
- Use the Quickclot Bleeding Control Dressing only as directed. Ensure that the area is clean before application.
- When using the Mini Compression Dressing, apply firm pressure to the wound to control bleeding.
- The Hypothermia Blanket should be used to maintain body temperature in cases of shock or exposure to cold.
- Flat Fold Tape should be used to secure dressings in place. Ensure it does not come into contact with open wounds.
- If you are unsure about how to use any of the supplies, seek guidance from a qualified medical professional.

### Instructions for Installation and Usage

#### 1. Accessing the Kit:

- The Micro TKN can be mounted on a belt or MOLLE system.
- To access the supplies, pull the BLIP featured pull tabs with one hand or a single finger from either side.

#### 2. Using Medical Supplies:

- Trauma Gloves: Put on gloves before touching any medical supplies.
- Quickclot Dressing: Open the packaging, apply to the wound, and hold in place for the recommended time.
- Mini Compression Dressing: Remove from packaging, place over the wound, and apply pressure.
- Hypothermia Blanket: Unfold the blanket and wrap it around the individual in need.
- Flat Fold Tape: Cut to the desired length and apply to secure dressings.

#### 3. PostUse:

- After using any supplies, dispose of them according to local regulations for medical waste.
- Restock the kit as necessary to ensure it is always ready for use.

## **Disposal Instructions**

- Dispose of used medical supplies in accordance with local regulations for hazardous waste.
- Ensure that the kit itself is disposed of properly if it is damaged or no longer needed.
- Recycle any components that are recyclable according to local guidelines.

## **Contact Information for Further Support**

For any inquiries regarding the Micro Trauma Kit NOW! Essentials Supplies, please refer to the product packaging for the relevant contact details. Ensure you have the product information ready when reaching out for support.

Thank you for prioritizing safety and preparedness with the Micro Trauma Kit NOW! Essentials Supplies.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)