

# MICRO TRAUMA KIT NOW! PRO SUPPLIES - BLUE FORCE GEAR

## MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT ODG

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW!, or Micro TKN, is the smallest version of the Trauma Kit NOW! for every day carry. It is designed to hold essential lifesaving medical supplies on a belt. The Micro TKN was designed as an Every Day Carry trauma kit for law enforcement professionals, prepared citizens, or hunters. Deployment of critical first aid supplies can be done with one hand or a single finger from either the left or right side by pulling the BLIP featured pull tabs. The Micro TKN consists of two main components - the outer MOLLE or Belt mounted pouch utilizing the Ten-Speed technology, and a removable insert that keeps medical supplies organized. The PRO Kit Includes: Trauma Gloves (1x Pair) Quickclot Combat Gauze LE (1x) SWAT-T Pressure Device (1x) Hypothermia Blanket - Green/Silver (1x)



### Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! PRO SUPPLIES MOLLE MOUNT ODG](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430100745
- Mfr. No.: HW-TKN-MTKN-PRO-OD
- Size: 1 Person
- Delivery weight: 0.222kg
- Shipping height: 51mm
- Shipping width: 152mm
- Shipping length: 229mm
- UPC: 810073650084

### Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [MICRO TRAUMA KIT NOW! PRO SUPPLIES Safety Instructions](#)
- [About Us](#)

# MICRO TRAUMA KIT NOW! PRO SUPPLIES Safety Instructions

## Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PRO SUPPLIES. This kit is designed to provide essential lifesaving medical supplies for everyday carry, ensuring you are prepared for emergencies. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

## General Safety Guidelines

- Ensure the Micro Trauma Kit NOW! is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure all supplies are present and within their expiration dates.
- This product is intended for adults and should not be used by children without adult supervision.
- Always follow the instructions provided with each medical supply included in the kit.
- In case of an emergency, seek professional medical help immediately, even if you administer first aid.

## Specific Safety Precautions for Use

- **Trauma Gloves:** Use gloves to protect yourself and the victim from infection. Dispose of gloves properly after use.
- **Quickclot Combat Gauze:** Follow the instructions for application carefully to ensure effective use. Do not apply to open wounds without proper training.
- **SWATT Pressure Device:** Ensure you understand how to use the pressure device correctly. Improper use can lead to further injury.
- **Hypothermia Blanket:** Use the blanket to prevent heat loss in victims. Ensure the victim is in a safe position before applying the blanket.
- **General Precautions:** Always wash your hands before and after using any medical supplies. If you are unsure about how to use a specific item, seek guidance from a qualified professional.

## Instructions for Installation and Usage

### • Accessing the Kit:

- The Micro TKN can be mounted on a belt or MOLLE system.
- To access supplies, pull the BLIP featured pull tabs using one hand or a single finger from either side.

### • Using the Medical Supplies:

- **Trauma Gloves:** Put on the gloves before administering first aid.
- **Quickclot Combat Gauze:**
  1. Open the package carefully.
  2. Apply the gauze to the wound with direct pressure.
  3. Secure it in place as needed.
- **SWATT Pressure Device:**
  1. Open the device and follow the instructions for applying pressure to the wound.
  2. Ensure the device is secured properly to maintain pressure.
- **Hypothermia Blanket:**
  1. Unfold the blanket and wrap it around the victim.
  2. Ensure that the shiny side is facing inwards to retain body heat.

- **After Use:**

- Dispose of any used medical supplies in accordance with local regulations.
- Replenish any used items in the kit to ensure it remains ready for future use.

## **Disposal Instructions**

- Dispose of used medical supplies according to your local waste management regulations.
- Do not dispose of medical supplies in regular household waste unless specified.
- Store the Micro Trauma Kit NOW! in a safe place until needed, and ensure it is kept out of reach of children.

## **Contact Information for Further Support**

For any inquiries or concerns regarding the MICRO TRAUMA KIT NOW! PRO SUPPLIES, please refer to the product packaging for contact details of the manufacturer or authorized distributor.

### **Important Note**

- In case of an emergency or if you encounter an unsafe product, report it to the relevant authorities.
- Stay informed on product recalls and safety updates through the EU's Safety Gate platform.

Thank you for prioritizing safety and preparedness with the MICRO TRAUMA KIT NOW! PRO SUPPLIES.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)