# TRAUMA KIT NOW! PRO SUPPLIES OD GREEN

The Trauma Kit NOW! – MEDIUM is a MOLLE mounted medical kit that can be deployed in seconds. This first aid pouch was originally designed to hold the contents of the Army IFAK and other medical supplies. This medical kit mounts on a 3x3 MOLLE field, a battle belt or other MOLLE Compatible Surfaces using the patented BFG Helium Whisper Attachment System. The Helium Whisper Attachment System is the lightest MOLLE attachment system in existence. Featuring less bulk so it sits closer to the body armor, pack or belt system. It has fewer parts to fail and is metal free to reduce secondary fragmentation. Best of all, Helium Whisper offers overall weight savings of up to 40% over traditional pouches. PRO SUPPLIES INCLUDE: Trauma Gloves (Rolled 1x Pair) Quikclot Bleeding Control Dressing (1x) Emergency Trauma Dressing 6" (1x) Hypothermia Blanket - Green / Silver (1x) Flat Fold Tape (1x) Trauma Shears (1x) Full Size Chest Seal (1x Twin Pack)

## Attributes

- Name: TRAUMA KIT NOW! PRO SUPPLIES OD GREEN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100750
- Mfr. No.: HW-TKN-MED-PRO-
- Size: 1 Person
- Delivery weight: 0.492kg
- Shipping height: 83mm
- Shipping width: 133mm
- Shipping length: 191mm
- UPC: 810073650282

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# Trauma Kit NOW! – Medium Safety Instruction Guide

## Introduction

Thank you for choosing the Trauma Kit NOW! – Medium. This kit is designed to provide essential medical supplies for emergency situations. It is important to read and understand this safety instruction guide to ensure safe and effective use of the product.

## General Safety Guidelines

- **Product Safety:** Ensure that the kit is used only for its intended purpose. This product meets the safety requirements outlined in the EU General Product Safety Regulation (GPSR).
- **Enhanced Recalls:** Stay informed about any recalls or safety notices related to this product. You can check for updates on the EU Safety Gate platform.
- **Online Shopping:** If you purchased this product online, ensure that the seller complies with safety regulations.
- **Special Consumer Focus:** Keep this kit out of reach of children and vulnerable individuals. Adult supervision is recommended when using the kit.
- **EU Contact Point:** For any safety inquiries or concerns, refer to the provided contact information in the product packaging.
- **Rapid Alerts:** Be aware that safety alerts regarding this product may be issued. Regularly check the Safety Gate system for updates.

## Specific Safety Precautions for Use

- Always inspect the kit before use to ensure all components are intact and undamaged.
- Use the kit only for its intended purpose. Do not use any items that are expired or damaged.
- If you are unsure about how to use any of the medical supplies, seek guidance from a qualified medical professional.
- Avoid using the kit in environments that may compromise its integrity, such as extreme temperatures or wet conditions.
- Ensure that the kit is easily accessible in case of an emergency.

## Instructions for Installation and Usage

### 1. Mounting the Kit:

- Ensure that the mounting surface is clean and dry.
- Use the MOLLE attachment system to securely fasten the kit to a 3x3 MOLLE field, battle belt, or compatible surface.
- Verify that the kit is securely attached before use.

### 2. Using the Medical Supplies:

- **Trauma Gloves:** Always wear gloves when providing first aid to prevent contamination.
- **Quikclot Bleeding Control Dressing:** Apply directly to the wound and apply pressure as needed.
- **Emergency Trauma Dressing:** Use this dressing for larger wounds and secure it in place.
- **Hypothermia Blanket:** Use to maintain body temperature in case of shock or exposure to cold.
- **Flat Fold Tape:** Use to secure dressings or bandages in place.
- **Trauma Shears:** Use to cut clothing or bandages as necessary.
- **Full Size Chest Seal:** Apply to any open chest wounds to prevent air from entering the chest cavity.

### 3. PostUse Care:

- Dispose of any used medical supplies according to local regulations.
- Replenish any used items in the kit to ensure it is ready for future use.

## Disposal Instructions

- Dispose of used medical supplies in accordance with local waste disposal regulations.
- Ensure that all items are disposed of in a manner that prevents harm to others or the environment.
- Recycle packaging materials where possible.

## **Contact Information for Further Support**

For any questions or concerns regarding the Trauma Kit NOW! – Medium, please refer to the contact information provided in the product packaging. It is important to have access to support for any safety inquiries or productrelated issues.

Thank you for choosing the Trauma Kit NOW! – Medium. Your safety is our priority.

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