MICRO TRAUMA KIT NOW! - NANO - MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-COYOTE BROWN

The Micro Trauma Kit NOW! Nano is the smallest of the TKN medical kits. This mini trauma pouch holds the essential medical gear intended to stop the number one cause of traumatic death: bleeding! This first aid pouch can easily be stowed in a pocket, medical pack, and other similar sized spaces. ESSENTIAL INCLUDES: Trauma Gloves (Rolled 1x Pair), Quikclot Bleeding Control Dressing (1x), and a Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), and a Pressure Device (1x) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), and a Pressure Device (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW!-NANO-ESSENTIAL SUPPLIES-COYOTE BROWN
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100987
- Mfr. No.: TKN-NANO-ESS-CB
- Size: 1 Person
- Delivery weight: 0.147kg
- Shipping height: 32mm
- Shipping width: 89mm
- Shipping length: 171mm
- UPC: 810073652828

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MICRO TRAUMA KIT NOW! NANO MICRO TRAUMA KIT NOW! Safety Instruction Guide

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! Nano. This compact trauma kit is designed to provide essential medical supplies to help manage traumatic injuries, particularly those involving bleeding. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- Always ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace any used or expired items promptly.
- Keep the kit out of reach of children and vulnerable individuals unless under adult supervision.
- In case of an emergency, call for professional medical help immediately after using the kit.
- Do not use the kit for purposes other than its intended use.
- Report any unsafe products or accidents to local authorities.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Always wear gloves when handling blood or bodily fluids to reduce the risk of infection.
 - Dispose of gloves properly after use to prevent contamination.
- Quikclot Bleeding Control Dressing:
 - Apply the dressing directly to the wound and apply pressure to control bleeding.
 - Do not remove the dressing once applied; it may interfere with clotting.
- Pressure Device:
 - Ensure the device is used as instructed to effectively control bleeding.
 - Do not apply excessive pressure, as this may cause further injury.

Instructions for Installation and Usage

- **Storage:**
 - Store the Micro Trauma Kit in a designated area that is easily accessible in case of emergencies.
- **Usage Steps:**
 1. Assess the situation and ensure safety for yourself and the injured person.
 2. Put on trauma gloves before attending to the injured person.
 3. If bleeding is present, use the Quikclot Bleeding Control Dressing:
 - Remove the dressing from its packaging.
 - Apply the dressing directly to the wound.
 - Use your hand to apply firm pressure for at least 3 minutes.
 4. If bleeding continues, apply a second dressing without removing the first.
 5. Use the Pressure Device if necessary, following the specific instructions provided with the device.
 6. Once bleeding is controlled, monitor the injured person and seek professional medical assistance.

Disposal Instructions

- Dispose of used gloves and any contaminated materials in accordance with local regulations for hazardous waste.
- Check local guidelines for disposal of medical supplies and follow them accordingly.

Contact Information for Further Support

For any inquiries or additional support, please reach out to the designated contact point for your region. It is important to have this information readily available for any safety concerns or product-related questions.

By following these guidelines, you can ensure the safe and effective use of your MICRO TRAUMA KIT NOW! Nano. Always prioritize safety and seek professional medical assistance when necessary. Thank you for your attention to these safety instructions.

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