MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ESSENTIAL -WOLF GRAY

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. **ESSENTIALS INCLUDES:** Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) **PRO INCLUDES:** Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) **ADVANCED INCLUDES:** Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)

Attributes

- Name: MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ESSENTIAL -WOLF GRAY
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430101005
- Mfr. No.: HW-TKN-PLUS-ESS
- Size: 1 Person
- Delivery weight: 0.458kg
- Shipping height: 76mm
- Shipping width: 203mm
- Shipping length: 254mm
- UPC: 810073653054

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MICRO TRAUMA KIT NOW! PLUS+ User Safety Instruction Guide

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+. This guide provides essential safety instructions to ensure the safe use and maintenance of your trauma kit. It is designed to help you understand the product's features, safety precautions, and proper usage in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and signs of damage.
- Keep the kit out of reach of children unless under adult supervision.
- Familiarize yourself with the contents and their intended use before an emergency occurs.
- Report any unsafe products or incidents related to this kit to the appropriate authorities.
- Stay informed about product recalls through the EU's Safety Gate platform.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling blood or bodily fluids to minimize infection risk.
- Use QuikClot products according to the provided instructions to effectively control bleeding.
- Ensure the emergency trauma dressing is applied firmly but not too tightly to avoid cutting off circulation.
- The nasopharyngeal airway should only be used by individuals trained in its application.
- The decompression needle is a medical device; use it only if trained to do so.
- If you are unsure about using any item in the kit, seek professional medical assistance immediately.

Instructions for Installation and Usage

Installation

- The kit can be mounted on MOLLE-compatible surfaces such as plate carriers, vests, or battle belts.
- To attach the kit, align the MOLLE straps with the webbing on your gear.
- Thread the straps through the webbing and secure them tightly to prevent movement.

Usage

1. **Assess the Situation:** Before using the kit, evaluate the scene for safety and the nature of injuries.
2. **Don Trauma Gloves:** Always wear the trauma gloves provided before administering first aid.
3. **Control Bleeding:**
 - Use the QuikClot Bleeding Control Dressing by applying it directly to the wound and applying firm pressure.
 - If bleeding does not stop, apply a second dressing and seek medical help.
4. **Apply Emergency Trauma Dressing:**
 - Place the dressing over the wound and secure it in place.
 - Monitor the patient for signs of shock and seek emergency medical assistance.
5. **Maintain Body Temperature:**
 - Use the hypothermia blanket to keep the patient warm and prevent shock.
6. **Airway Management:**
 - If the patient is unconscious and not breathing, consider using the nasopharyngeal airway as per your training.
7. **Decompression:**
 - If necessary and trained, use the decompression needle for tension pneumothorax.
8. **Chest Seal Application:**
 - Apply the fullsize chest seal over open chest wounds to prevent air from entering the chest cavity.

Disposal Instructions

- Dispose of used or expired items in accordance with local regulations for medical waste.
- Ensure that all sharp objects, such as needles, are disposed of in a designated sharps container.
- Do not dispose of the kit in regular household waste.

Contact Information for Further Support

For any questions or concerns regarding the MICRO TRAUMA KIT NOW! PLUS+, please consult the manufacturer's website or customer service for assistance. Ensure you have the product details and any relevant information ready when seeking support.

By following these safety instructions, you can ensure the effective use of your MICRO TRAUMA KIT NOW! PLUS+ while adhering to safety regulations. Always prioritize safety and seek professional help when necessary.

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