

MICRO TRAUMA KIT NOW! - PLUS+ - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO - BLACK

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO - BLACK](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101006
- Mfr. No.: HW-TKN-PLUS-PRO-BK
- Size: 1 Person
- Delivery weight: 0.245kg
- Shipping height: 102mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653061

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [MICRO TRAUMA KIT NOW! PLUS+ Safety Instruction Guide](#)
- [About Us](#)

MICRO TRAUMA KIT NOW! PLUS+ Safety Instruction Guide

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+ by BLUE FORCE GEAR. This trauma kit is designed to provide essential medical supplies for emergency situations. This guide offers important safety instructions to ensure the safe use and maintenance of the kit, in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Ensure that all users of the MICRO TRAUMA KIT NOW! PLUS+ are trained in first aid and emergency response.
- Regularly check the contents of the kit to ensure that all items are present and in good condition.
- Store the kit in a cool, dry place away from direct sunlight to maintain the integrity of the medical supplies.
- Keep the kit out of reach of children unless under adult supervision.
- Dispose of any expired or damaged items according to local regulations.

Specific Safety Precautions for Use

- Always wear trauma gloves when handling blood or bodily fluids to protect yourself and others.
- Use QuikClot products only as directed to avoid complications.
- Ensure that the emergency trauma dressing is applied correctly to prevent further injury.
- Be cautious when using sharp instruments, such as the decompression needle, to avoid accidental injury.
- If you are unsure about how to use any item in the kit, seek professional guidance before attempting to provide care.

Instructions for Installation and Usage

• Mounting the Kit:

- Use the MOLLE med kit mount to secure the kit to your plate carrier, vest, or battle belt.
- Ensure that the kit is firmly attached and easily accessible in case of an emergency.

• Using the Kit:

- In case of a traumatic injury, assess the situation and ensure your safety first.
- Put on trauma gloves before providing care.
- Follow these steps for key items:
 - **QuikClot Bleeding Control Dressing:**
 - Open the packaging carefully.
 - Apply directly to the wound, applying pressure as necessary.
 - **Emergency Trauma Dressing:**
 - Open the dressing and place it over the wound.
 - Secure it in place using the attached ties.
 - **Chest Seal:**
 - Remove the chest seal from its packaging.
 - Apply it over any open chest wounds to prevent air from entering the chest cavity.
 - **Hypothermia Blanket:**
 - Unfold the blanket and wrap it around the injured person to maintain body temperature.

- **After Use:**

- Dispose of used items according to local regulations.
- Restock the kit as necessary to ensure it is always ready for use.

Disposal Instructions

- Dispose of expired or damaged items in accordance with local waste disposal regulations.
- For items that have come into contact with blood or bodily fluids, follow biohazard disposal guidelines.
- Recycle packaging materials where possible.

Contact Information for Further Support

For any inquiries regarding the MICRO TRAUMA KIT NOW! PLUS+ or to report unsafe products, please refer to your local consumer safety authority or visit the EU's Safety Gate platform for updates on recalls and safety alerts.

Thank you for your attention to these safety guidelines. Your safety and the safety of others depend on the proper use and maintenance of the MICRO TRAUMA KIT NOW! PLUS+.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk