

MICRO TRAUMA KIT NOW! - PLUS+ - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO -RANGER GREEN

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO -RANGER GREEN](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101009
- Mfr. No.: HW-TKN-PLUS-PRO-RG
- Size: 1 Person
- Delivery weight: 0.236kg
- Shipping height: 305mm
- Shipping width: 83mm
- Shipping length: 254mm
- UPC: 810073653092

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [MICRO TRAUMA KIT NOW! PLUS+ Safety Instruction Guide](#)
- [About Us](#)

MICRO TRAUMA KIT NOW! PLUS+ Safety Instruction Guide

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+ by BLUE FORCE GEAR. This medical pouch is designed for individuals who require quick access to essential medical supplies in emergency situations. It is crucial to use this product safely and effectively to ensure your wellbeing and that of others.

General Safety Guidelines

- Always read and understand the instructions before using the product.
- Keep the kit out of reach of children unless they are under adult supervision.
- Regularly check the contents of the kit to ensure all items are present and in good condition.
- Do not use expired medical supplies.
- If you experience any adverse reactions or complications while using the kit, seek medical assistance immediately.
- Report any unsafe products or accidents to the relevant authorities.

Specific Safety Precautions for Use

- Ensure that you are trained in basic first aid before using the contents of this kit.
- Use trauma gloves to prevent contamination and protect yourself when treating injuries.
- Be cautious when using sharp objects such as scissors or needles included in the kit.
- Ensure that all dressings and gauze are applied according to the instructions provided with each item.
- For individuals with allergies, check the packaging of each item for potential allergens.
- If you are unsure about how to use any item in the kit, consult a medical professional or refer to first aid resources.

Instructions for Installation and Usage

1. Mounting the Kit:

- Attach the MOLLE med kit mount securely to your plate carrier, vest, or battle belt.
- Ensure that the kit is easily accessible for quick retrieval during emergencies.

2. Using the Kit:

- In the event of an emergency, calmly assess the situation and the injuries.
- Put on trauma gloves before handling any medical supplies.
- Use the QuikClot Bleeding Control Dressing to manage severe bleeding.
- Apply the Emergency Trauma Dressing as needed, ensuring it is secured properly.
- Utilize the Flat Fold Tape to secure dressings in place.
- If necessary, use the Full Size Chest Seal to treat chest wounds.

3. After Use:

- Dispose of any used medical supplies according to local regulations.
- Replenish any items that have been used from the kit to maintain readiness for future emergencies.
- Store the kit in a cool, dry place away from direct sunlight.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations for hazardous waste.
- Ensure that all items are securely packaged before disposal to prevent any risk of contamination.
- Check local guidelines for proper disposal of medical waste to ensure safety for the community.

Contact Information for Further Support

For any inquiries regarding safety or product information, please refer to the EU contact point or consult the product's official website. Always ensure you have the latest information regarding recalls or safety updates.

By following these safety instructions, you can help ensure the effective and safe use of the MICRO TRAUMA KIT NOW! PLUS+. Your safety and the safety of others depend on proper preparation and knowledge. Thank you for prioritizing safety.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk