

# MICRO TRAUMA KIT NOW! - PLUS+ - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO -WOLF GRAY

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



## Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - PRO -WOLF GRAY](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101010
- Mfr. No.: HWTKNPLUSPRO
- Size: 1 Person
- Delivery weight: 0.503kg
- Shipping height: 51mm
- Shipping width: 203mm
- Shipping length: 254mm
- UPC: 810073653108

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [MICRO TRAUMA KIT NOW! PLUS+ User Safety Instruction Guide](#)
- [About Us](#)

# MICRO TRAUMA KIT NOW! PLUS+ User Safety Instruction Guide

## Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+ from BLUE FORCE GEAR. This guide provides essential safety instructions and information to ensure the safe and effective use of your trauma kit. Please read this guide carefully and keep it for future reference.

## General Safety Guidelines

- Ensure that the kit is used only for its intended purpose as a first aid and trauma response tool.
- Always check the contents of the kit regularly to ensure that all items are present and within their expiration dates.
- Store the kit in a cool, dry place, away from direct sunlight, to maintain the integrity of the medical supplies.
- Be aware of the potential hazards associated with the use of medical supplies, including allergic reactions and the risk of infection.
- Keep the kit out of reach of children and vulnerable individuals unless supervised by an adult.
- Report any unsafe products or incidents related to the kit to the relevant authorities.

## Specific Safety Precautions for Use

- Always wear gloves when administering first aid to prevent crosscontamination and protect yourself and the patient.
- Ensure that all items, especially dressings and gauze, are sterile before use.
- Follow the instructions provided for each item carefully to avoid misuse and ensure effective treatment.
- If you are unsure about how to use any item in the kit, seek assistance from a trained medical professional.
- Be cautious when using sharp objects, such as scissors and needles, to avoid injury.
- Keep the kit away from sources of heat and moisture to prevent damage to the contents.

## Instructions for Installation and Usage

### 1. Familiarize Yourself with the Kit:

- Review the contents of the MICRO TRAUMA KIT NOW! PLUS+ to understand what items are included and their specific uses.

### 2. Using the Kit:

- In case of an emergency, assess the situation and ensure it is safe to proceed.
- Use the trauma gloves provided to protect yourself and the patient.
- Administer first aid using the appropriate items from the kit based on the injury type.
- For bleeding control, use the QuikClot Bleeding Control Dressing as per the instructions on the packaging.
- If a chest seal is needed, follow the instructions for applying the Full Size Chest Seal carefully.
- In case of a severe injury, do not hesitate to call emergency services for professional assistance.

### 3. PostUsage Care:

- Dispose of used items, especially gloves and dressings, in a safe manner to prevent contamination.
- Restock any used items promptly to ensure the kit is ready for future emergencies.

## **Disposal Instructions**

- Dispose of any expired or used medical supplies in accordance with local regulations for medical waste.
- Do not throw away medical items in regular household waste unless specified by local guidelines.
- Ensure that any sharp objects, such as needles, are disposed of in a designated sharps container.

## **Contact Information for Further Support**

For safety inquiries or assistance regarding the MICRO TRAUMA KIT NOW! PLUS+, please refer to the contact information provided on the packaging or the manufacturer's website.

## **Conclusion**

The MICRO TRAUMA KIT NOW! PLUS+ is designed to provide essential first aid in emergencies. By following the safety instructions outlined in this guide, you can ensure the safe and effective use of your trauma kit. Stay prepared and prioritize safety at all times.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)