

# MICRO TRAUMA KIT NOW! - PLUS+ - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ADVANCED - MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



## Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ADVANCED - MULTICAM](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101013
- Mfr. No.: HW-TKN-PLUS-ADV-MC
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 305mm
- Shipping width: 83mm
- Shipping length: 254mm
- UPC: 810073653139

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [Micro Trauma Kit NOW! PLUS+ Safety Instructions](#)
- [About Us](#)

# Micro Trauma Kit NOW! PLUS+ Safety Instructions

## Introduction

Thank you for choosing the Micro Trauma Kit NOW! PLUS+ by Blue Force Gear. This kit is designed for quick and effective trauma response. It is essential to follow the safety guidelines and instructions provided in this manual to ensure safe use and optimal performance.

## General Safety Guidelines

- **Product Safety:** This product is designed to be safe for all users when used as intended. Ensure that you are familiar with its contents and their purposes.
- **Enhanced Recalls:** Stay informed about any product recalls or safety notices. If you receive a recall notice, follow the instructions provided.
- **Online Shopping:** If purchased online, ensure that the retailer complies with safety regulations.
- **Special Consumer Focus:** This product is not intended for use by children. Keep it out of reach of minors.
- **EU Contact Point:** For any safety inquiries, please refer to the contact information provided by your retailer.
- **Rapid Alerts:** Monitor the EU's Safety Gate platform for any updates on product safety.

## Specific Safety Precautions for Use

- **Read Instructions:** Always read the instructions carefully before using the kit.
- **Check Contents:** Regularly check the contents of the kit to ensure all items are present and in good condition.
- **Use as Intended:** This kit is designed for trauma situations. Do not use it for purposes other than first aid.
- **Personal Protective Equipment:** Always wear trauma gloves when administering first aid to prevent infection.
- **Avoid Contamination:** Ensure that all medical supplies are used in a clean environment to avoid contamination.
- **Seek Professional Help:** In the event of a serious injury, always seek professional medical assistance as soon as possible.

## Instructions for Installation and Usage

### 1. Accessing the Kit:

- The kit should be easily accessible in emergencies. Store it in a location known to all potential users.
- Ensure that the MOLLE med kit mount is securely attached to the plate carrier, vest, or battle belt.

### 2. Using the Contents:

- **Trauma Gloves:** Roll out a pair of gloves before touching any injured person.
- **QuikClot Bleeding Control Dressing:** Apply directly to the wound and secure it in place.
- **Emergency Trauma Dressing:** Use the dressing to cover larger wounds. Follow the instructions on the packaging for proper application.
- **Hypothermia Blanket:** Use to keep the injured person warm.
- **Flat Fold Tape:** Use for securing dressings and other medical supplies.
- **Pressure Device:** Apply to control severe bleeding.

- **Chest Seal:** Use as instructed to treat penetrating chest wounds.
- **Nasopharyngeal Airway:** Follow the instructions for correct placement.
- **Decompression Needle:** Only use if trained; improper use can cause harm.

### 3. **PostUse:**

- Dispose of used medical supplies properly according to local regulations.
- Restock the kit as necessary to ensure it is always ready for use.

## **Disposal Instructions**

- Dispose of all used medical supplies in accordance with local waste disposal regulations.
- Do not dispose of medical items in regular household waste.
- Ensure that all sharp objects, such as needles, are disposed of in a designated sharps container.

## **Contact Information for Further Support**

For any questions or concerns regarding the Micro Trauma Kit NOW! PLUS+, please refer to the contact information provided by your retailer or manufacturer. They can assist with safety inquiries and provide additional support as needed.

By following these safety instructions, you can ensure the effective use of the Micro Trauma Kit NOW! PLUS+ while maintaining safety for yourself and others. Thank you for prioritizing safety and preparedness.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)