

MICRO TRAUMA KIT NOW! - PLUS+ - BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ADVANCED -WOLF GRAY

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Micro Trauma Kit NOW! - PLUS+ is 30% larger than the original BFG micro medic kit (MTKN). This medical pouch is for those users who want the design of the MTKN but need to carry a full-size pressure dressing and chest seal. Available with a MOLLE med kit mount that adapts to plate carriers, vests, battle belts or other MOLLE compatible surfaces. It offers the same high speed access to your medical gear as the other BFG MED MTKN and TKN trauma first aid pouches. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Bleeding Control Dressing (1x), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), and Pressure Device (1x) PRO INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Pressure Device (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 1x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR MICRO TRAUMA KIT NOW! - PLUS+ - MOLLE - ADVANCED -WOLF GRAY](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101015
- Mfr. No.: HW-TKN-PLUS-ADV-WF
- Size: 1 Person
- Delivery weight: 0.503kg
- Shipping height: 51mm
- Shipping width: 203mm
- Shipping length: 254mm
- UPC: 810073653153

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [MICRO TRAUMA KIT NOW! PLUS+ Safety Instruction Guide](#)
- [About Us](#)

MICRO TRAUMA KIT NOW! PLUS+ Safety Instruction Guide

Introduction

Thank you for choosing the MICRO TRAUMA KIT NOW! PLUS+ by BLUE FORCE GEAR. This advanced medical pouch is designed to provide essential first aid supplies in a compact and accessible format. It is important to follow the safety guidelines outlined in this document to ensure the safe and effective use of this product.

General Safety Guidelines

- Ensure that all users are trained in first aid and emergency response before using the Micro Trauma Kit NOW! PLUS+.
- Regularly check the contents of the kit to ensure that all items are present, undamaged, and within their expiration dates.
- Keep the kit out of reach of children unless under adult supervision.
- Store the kit in a cool, dry place, away from direct sunlight and extreme temperatures.
- Report any unsafe products or incidents to the appropriate authorities.

Specific Safety Precautions for Use

- Always wear trauma gloves when administering first aid to prevent contamination and infection.
- Be aware of the potential for allergic reactions to any materials in the kit, such as latex gloves or adhesive dressings.
- Use the QuikClot products according to the instructions provided, ensuring that they are applied correctly to control bleeding.
- Handle the decompression needle with care to prevent injury. Ensure that you are trained in its use.
- If using the nasopharyngeal airway, ensure proper sizing and insertion technique to avoid injury.

Instructions for Installation and Usage

1. Mounting the Kit:

- Attach the Micro Trauma Kit NOW! PLUS+ to MOLLEcompatible surfaces, such as plate carriers or battle belts, using the provided mounting system.
- Ensure that the kit is securely fastened and easily accessible.

2. Using the Kit:

- Open the kit and familiarize yourself with the contents.
- In case of an emergency, quickly locate the necessary supplies.
- Follow these steps for specific items:
 - **Trauma Gloves:** Roll and wear them before providing first aid.
 - **QuikClot Bleeding Control Dressing:** Apply directly to the wound and hold in place to control bleeding.
 - **Emergency Trauma Dressing:** Place over the wound and secure with the provided tape.
 - **Hypothermia Blanket:** Wrap the victim to prevent heat loss.
 - **Flat Fold Tape:** Use to secure dressings or splints.
 - **Pressure Device:** Apply to control severe bleeding.
 - **Chest Seal:** Follow instructions for application to treat chest wounds.
 - **Nasopharyngeal Airway:** Ensure proper sizing and insert gently.
 - **Decompression Needle:** Use only if trained; follow emergency protocols.

3. **PostUse:**

- Dispose of used items according to local regulations.
- Replenish any used supplies to maintain the kit's readiness for future emergencies.

Disposal Instructions

- Dispose of expired or damaged items in accordance with local waste disposal regulations.
- Used medical supplies, such as gloves and dressings, should be placed in biohazard waste containers where applicable.
- Recycle packaging materials where possible.

Contact Information for Further Support

For any inquiries regarding the MICRO TRAUMA KIT NOW! PLUS+, including safety concerns and product support, please refer to your local distributor or visit the official BLUE FORCE GEAR website for assistance.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk