

# TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES - COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



## Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES - COYOTE](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101017
- Mfr. No.: HW-TKN-SMALL-ES
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 76mm
- Shipping width: 229mm
- Shipping length: 254mm
- UPC: 810073653221

## Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [Trauma Kit NOW! SMALL Safety Instructions](#)
- [About Us](#)

# Trauma Kit NOW! SMALL Safety Instructions

## Introduction

Thank you for choosing the Trauma Kit NOW! SMALL by Blue Force Gear. This medical kit is designed to provide essential supplies for emergency situations. To ensure safe and effective use, please read and follow these safety instructions carefully.

## General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace any expired items.
- Keep the kit out of reach of children and vulnerable individuals unless supervised by an adult.
- Use the kit only for its intended purpose—emergency first aid.
- In case of a serious medical emergency, seek professional medical assistance immediately.

## Specific Safety Precautions for Use

- Always wear the Trauma Gloves provided when handling blood or bodily fluids to reduce the risk of contamination.
- Ensure that the QuikClot and other bleeding control products are applied according to the instructions provided in the kit.
- Be cautious when using sharp items, such as scissors or needles, to avoid injury.
- If using the Nasopharyngeal Airway, ensure it is the correct size for the individual and follow the insertion instructions carefully.
- Do not attempt to use the Decompression Needle unless you are trained in its application.

## Instructions for Installation and Usage

### • Opening the Kit:

- Open the MOLLE pouch by unfastening the straps.
- Ensure all items are accessible and in good condition.

### • Using Trauma Gloves:

- Pull out the rolled gloves and wear them before handling any supplies.
- Dispose of gloves properly after use to maintain hygiene.

### • Applying QuikClot:

- Open the package and apply directly to the bleeding wound.
- Apply pressure to the dressing for at least 3 minutes or until bleeding stops.

### • Using the Emergency Trauma Dressing:

- Remove the dressing from its packaging.
- Place the dressing over the wound and secure it in place using the attached ties.

### • Using the Hypothermia Blanket:

- Unfold the blanket and wrap it around the individual to maintain body heat.

### • Using Compressed Gauze:

- Open the gauze package and apply it to the wound, pressing firmly to control bleeding.
- **Using the Nasopharyngeal Airway:**
  - Lubricate the airway before insertion.
  - Gently insert the airway into the nostril, aiming towards the back of the throat.
- **Using the Decompression Needle:**
  - Only use if trained. Locate the appropriate site on the chest and insert the needle as directed.

## **Disposal Instructions**

- Dispose of used gloves, dressings, and any other single-use items in accordance with local regulations for medical waste.
- Ensure that sharp items, such as needles, are disposed of in a designated sharps container.
- Recycle any packaging materials where possible.

## **Contact Information for Further Support**

For questions or concerns regarding the Trauma Kit NOW! SMALL, please refer to the manufacturer's website or customer support for additional information.

Please remember to keep this safety instruction guide with your Trauma Kit NOW! SMALL for easy reference. Your safety and wellbeing are our top priorities.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)