

TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES-MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - ESSENTIALS SUPPLIES-MULTICAM](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101018
- Mfr. No.: HW-TKN-SMALL-ES
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 279mm
- Shipping width: 229mm
- Shipping length: 76mm
- UPC: 810073653238

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [TRAUMA KIT NOW! SMALL BLUE FORCE GEAR TRAUMA KIT NOW! SMALL MOLLE ESSENTIALS SUPPLIESMULTICAM Safety Instruction Guide](#)
- [About Us](#)

TRAUMA KIT NOW! SMALL BLUE FORCE GEAR TRAUMA KIT NOW! SMALL MOLLE ESSENTIALS SUPPLIESMULTICAM Safety Instruction Guide

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL from Blue Force Gear. This medical kit is designed for effective emergency response and is suitable for individual use. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

General Safety Guidelines

- **Product Safety:** Ensure safe use of the Trauma Kit by following all instructions and guidelines. This kit is designed for nonfood use and addresses risks associated with emergency medical situations.
- **Enhanced Recalls:** Stay informed about potential recalls. If a safety issue arises, standardized recall notices will be provided, and remedies will be offered at no cost.
- **Online Shopping:** If purchased online, ensure that the platform complies with safety requirements similar to those of physical stores.
- **Special Consumer Focus:** Extra precautions should be taken when this kit is used by or around vulnerable groups, such as children.
- **EU Contact Point:** For safety inquiries, please refer to the contact information provided with the product packaging.
- **Rapid Alerts:** Be aware of updates regarding unsafe products through the EU's Safety Gate system.

Specific Safety Precautions for Use

- **Age Restrictions:** This kit is intended for adult use. Keep out of reach of children.
- **Check Contents:** Before use, ensure all components are present and in good condition. Do not use if any item is damaged or missing.
- **Training:** Familiarize yourself with basic first aid procedures. Proper training can enhance the effectiveness of the kit.
- **Hygiene:** Use gloves provided in the kit to maintain hygiene and prevent infection during treatment.
- **Storage:** Store the kit in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Expiry Dates:** Regularly check the expiry dates of the medical supplies included in the kit. Replace expired items promptly.

Instructions for Installation and Usage

1. Preparation:

- Open the pouch and lay out all contents to ensure you have everything needed for a response.
- Familiarize yourself with the use of each item included in the kit.

2. Using the Kit:

- **Trauma Gloves:** Always wear gloves when providing first aid to prevent contamination.
- **QuikClot Bleeding Control Dressing:** Apply directly to the wound and apply pressure to control bleeding.
- **Emergency Trauma Dressing:** Use as needed for severe bleeding. Follow the instructions on the packaging.
- **Hypothermia Blanket:** Use to maintain body temperature in case of shock or exposure.
- **Flat Fold Tape and Compressed Gauze:** Use for securing dressings and gauze to the wound.

site.

- **Chest Seal:** Apply over open chest wounds to prevent air from entering the chest cavity.
- **Nasopharyngeal Airway and Decompression Needle:** These items should only be used by trained individuals.

3. **PostUse:**

- Dispose of used items safely according to local regulations.
- Restock any used supplies as soon as possible to maintain the kit's readiness.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations.
- Do not dispose of items in regular household waste. Follow guidelines for hazardous waste if applicable.
- Ensure the kit is kept out of reach of children during disposal to prevent accidental use.

Contact Information for Further Support

For any questions or concerns regarding the use of the TRAUMA KIT NOW! SMALL, please refer to the contact information provided with your product packaging. It is important to reach out for assistance if needed to ensure safe and effective use of your trauma kit.

Thank you for your attention to these safety instructions. Stay safe and prepared!

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk