

TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL-MOLLE-ESSENTIALS SUPPLIES-RANGER GREEN

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL-MOLLE-ESSENTIALS SUPPLIES-RANGER GREEN](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101019
- Mfr. No.: HW-TKN-SMALL-ES
- Size: 1 Person
- Delivery weight: 0.454kg
- Shipping height: 51mm
- Shipping width: 229mm
- Shipping length: 305mm
- UPC: 810073653245

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [TRAUMA KIT NOW! SMALL User Safety Instruction Guide](#)
- [About Us](#)

TRAUMA KIT NOW! SMALL User Safety Instruction Guide

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL by Blue Force Gear. This compact trauma kit is designed to provide essential medical supplies in a convenient, easyto carry format. Please read this safety instruction guide carefully to ensure safe and effective use of your trauma kit.

General Safety Guidelines

- Ensure that the trauma kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure all items are present and in good condition.
- Keep the trauma kit out of reach of children to prevent accidental use.
- Familiarize yourself with the contents and their uses before an emergency arises.
- In case of an emergency, always call emergency services first before administering any aid.
- Follow the instructions provided for each item in the kit to avoid misuse or injury.

Specific Safety Precautions for Use

- **Trauma Gloves:** Wear gloves when providing first aid to prevent contamination and protect both the caregiver and the patient.
- **QuikClot Products:** Use as directed. Avoid contact with eyes and mucous membranes. If contact occurs, rinse thoroughly with water.
- **Emergency Trauma Dressing:** Ensure the dressing is applied firmly but not too tight to avoid restricting blood flow.
- **Hypothermia Blanket:** Use to maintain body temperature in a patient who is experiencing shock or hypothermia.
- **Compressed Gauze:** Apply directly to wounds to control bleeding. Do not remove if it becomes soaked; apply additional gauze on top if necessary.
- **Chest Seal:** Ensure proper placement over open chest wounds. Do not remove once applied unless instructed by medical professionals.

Instructions for Installation and Usage

1. Opening the Kit:

- Unzip or unbuckle the trauma kit carefully to avoid damaging the contents.
- Lay the kit flat on a stable surface for easy access.

2. Using the Supplies:

- **Trauma Gloves:** Roll out the gloves and wear them before touching the patient or any medical supplies.
- **QuikClot Bleeding Control Dressing:** Open the packaging and apply directly to the wound. Apply pressure as needed.
- **Emergency Trauma Dressing:** Open the package, remove the dressing, and apply it directly to the injury. Secure it in place.
- **Hypothermia Blanket:** Unfold and wrap the patient to retain heat. Ensure the shiny side is facing the patient for optimal insulation.
- **Flat Fold Tape:** Use to secure dressings in place as needed.
- **Compressed Gauze:** Open the package and apply to the wound, applying pressure to control bleeding.

- **Chest Seal:** Follow the instructions on the packaging for proper application over the wound.
- **Nasopharyngeal Airway and Decompression Needle:** These items should only be used by trained personnel.

3. **PostUse:**

- Dispose of used items according to local regulations.
- Replenish any used supplies promptly to ensure the kit is always ready for use.

Disposal Instructions

- Dispose of used medical supplies in accordance with local regulations for hazardous waste.
- Do not throw away items in regular trash unless specified as safe.
- Ensure that any sharp objects (e.g., used needles) are disposed of in a designated sharps container.

Contact Information for Further Support

For any inquiries regarding the TRAUMA KIT NOW! SMALL, including safety concerns or product recalls, please refer to the contact information provided with your purchase or consult the Blue Force Gear website.

Conclusion

The TRAUMA KIT NOW! SMALL is an essential tool for emergency preparedness. By following these safety guidelines, you can ensure the safe and effective use of your trauma kit. Always prioritize safety and seek professional medical assistance when needed. Thank you for choosing Blue Force Gear.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk