

TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL-MOLLE-ESSENTIALS SUPPLIES-WOLF GRAY

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL-MOLLE-ESSENTIALS SUPPLIES-WOLF GRAY](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101020
- Mfr. No.: HW-TKN-SMALL-ESS-WF
- Size: 1 Person
- Delivery weight: 0.442kg
- Shipping height: 102mm
- Shipping width: 178mm
- Shipping length: 254mm
- UPC: 810073653252

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Trauma Kit NOW! SMALL Safety Instruction Guide](#)
- [About Us](#)

Trauma Kit NOW! SMALL Safety Instruction Guide

Introduction

Thank you for choosing the Trauma Kit NOW! SMALL from Blue Force Gear. This compact and essential medical kit is designed to provide critical supplies for emergency situations. To ensure safe and effective use of this product, please read and follow the safety guidelines outlined in this manual.

General Safety Guidelines

- **Product Safety:** This kit is designed for safe use during emergency situations. Always ensure that you understand how to use each item properly before an emergency occurs.
- **Enhanced Recalls:** In the event of a product recall, you will receive standardized notices and remedies at no cost.
- **Online Shopping:** If purchased online, you have the same safety protections as if you bought it in a physical store.
- **Special Consumer Focus:** This kit is not intended for use by children. Keep out of reach of minors.
- **EU Contact Point:** For safety inquiries, please refer to the contact information provided with your product.
- **Rapid Alerts:** Stay informed about potential safety issues by checking the EU's Safety Gate platform.

Specific Safety Precautions for Use

- Always check the expiration dates on all medical supplies before use.
- Inspect the kit regularly to ensure all items are intact and properly stored.
- Do not use any items that appear damaged, discolored, or compromised.
- Ensure that you are trained in basic first aid and the use of the items within this kit.
- Keep the kit in a dry, cool place to maintain the integrity of the supplies.
- In case of severe injuries, seek professional medical assistance immediately.

Instructions for Installation and Usage

1. Unpacking the Kit:

- Carefully remove the Trauma Kit NOW! SMALL from its packaging.
- Ensure all components listed in the product description are present.

2. Understanding the Contents:

- Familiarize yourself with each item in the kit:
 - Trauma Gloves (Rolled 2x Pair)
 - QuikClot Bleeding Control Dressing (1x white for Essentials, black for PR, green for Advanced)
 - Emergency Trauma Dressing 6" (1x)
 - Hypothermia Blanket (1x)
 - Flat Fold Tape (1x)
 - Compressed Gauze (1x)
 - Full Size Chest Seal (1x Twin Pack for PR and Advanced)
 - Nasopharyngeal Airway (1x for Advanced)
 - Decompression Needle (1x for Advanced)

3. Using the Kit:

- In an emergency, follow these steps:
 - Put on the trauma gloves to protect yourself and the injured person.

- Assess the situation and determine the necessary supplies.
- Use the appropriate dressing for bleeding control (QuikClot or Emergency Trauma Dressing).
- Apply the hypothermia blanket if the person is in shock or has lost a significant amount of blood.
- Utilize the flat fold tape and compressed gauze as needed for wound care.
- For advanced care, use the chest seal, nasopharyngeal airway, and decompression needle as per your training.

4. PostUse Care:

- After use, dispose of any singleuse items according to local regulations.
- Replenish any used supplies to maintain the kit's readiness for future emergencies.

Disposal Instructions

- Dispose of used or expired items in accordance with local regulations.
- For items that are not biodegradable, please use appropriate waste disposal methods to minimize environmental impact.

Contact Information for Further Support

For safety inquiries or further assistance, please refer to the contact information provided with your product packaging or visit the manufacturer's website.

By following these guidelines, you can ensure the safe and effective use of your Trauma Kit NOW! SMALL. Always be prepared and prioritize safety in emergency situations.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk