

TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - PRO SUPPLIES - COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - PRO SUPPLIES - COYOTE](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101022
- Mfr. No.: HW-TKN-SMALL-PR
- Size: 1 Person
- Delivery weight: 0.544kg
- Shipping height: 76mm
- Shipping width: 229mm
- Shipping length: 254mm
- UPC: 810073653276

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [TRAUMA KIT NOW! SMALL USER SAFETY INSTRUCTIONS](#)
- [About Us](#)

TRAUMA KIT NOW! SMALL USER SAFETY INSTRUCTIONS

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL from Blue Force Gear. This compact medical pouch is designed to provide essential trauma supplies in an easyto carry format. It is crucial to understand how to use this kit safely and effectively to ensure the best outcomes in emergency situations.

General Safety Guidelines

- Always inspect the trauma kit before use to ensure all components are intact and within their expiration dates.
- Keep the trauma kit out of reach of children to prevent accidental use or ingestion of items.
- Store the kit in a cool, dry place to maintain the integrity of the medical supplies.
- Familiarize yourself with the contents of the kit and their intended uses before an emergency occurs.
- In case of an emergency, seek professional medical assistance as soon as possible, even after using the kit.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Ensure gloves are free from tears or punctures before use.
 - Dispose of gloves properly after use to prevent contamination.
- QuikClot Bleeding Control Dressing:
 - Apply directly to the wound and maintain pressure as instructed.
 - Do not remove the dressing once applied; seek medical help.
- Emergency Trauma Dressing:
 - Follow the instructions for application and secure tightly to control bleeding.
- Hypothermia Blanket:
 - Use to prevent heat loss in injured persons, especially in cold environments.
- Flat Fold Tape:
 - Ensure the skin is clean and dry before application to secure dressings.
- Compressed Gauze:
 - Use to pack wounds as needed; avoid excessive pressure that may cause further injury.

Instructions for Installation and Usage

1. Opening the Kit:

- Locate the pull tab on the top of the pouch and pull to access the contents.
- Ensure you have a clear space to work when accessing the kit.

2. Using Trauma Gloves:

- Carefully remove a pair of gloves from the kit.
- Put on gloves before handling any supplies or assisting an injured person.

3. Applying QuikClot Bleeding Control Dressing:

- Open the dressing package.
- Place the dressing directly over the wound and apply firm pressure for at least 3 minutes.

4. Applying Emergency Trauma Dressing:

- Open the dressing and place it over the wound.
- Secure it tightly using the builtin ties.

5. Using the Hypothermia Blanket:

- Unfold the blanket and wrap it around the injured person, ensuring full coverage.

6. Using Flat Fold Tape:

- Cut a piece of tape to the required length.
- Apply the tape to secure dressings or gauze in place.

7. Using Compressed Gauze:

- Remove the gauze from the kit.
- Pack the gauze into the wound as needed and secure with tape.

Disposal Instructions

- Dispose of used medical supplies, including gloves, dressings, and gauze, in accordance with local regulations.
- Do not reuse any disposable items; they are meant for single use only.
- Follow local guidelines for the disposal of hazardous waste if applicable.

Contact Information for Further Support

For safety inquiries or additional support, please refer to the contact information provided by the manufacturer. Ensure to have your product details ready for assistance.

By following these safety instructions, you can ensure the effective use of the TRAUMA KIT NOW! SMALL in emergencies while adhering to the EU General Product Safety Regulation (GPSR). Always prioritize safety and seek professional medical help when needed.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk