

TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL-MOLLE- PRO SUPPLIES-WOLF GRAY

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL-MOLLE- PRO SUPPLIES-WOLF GRAY](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101025
- Mfr. No.: HW-TKN-SMALL-PR
- Size: 1 Person
- Delivery weight: 0.442kg
- Shipping height: 102mm
- Shipping width: 178mm
- Shipping length: 254mm
- UPC: 810073653306

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [TRAUMA KIT NOW! SMALL USER SAFETY INSTRUCTIONS](#)
- [About Us](#)

TRAUMA KIT NOW! SMALL USER SAFETY INSTRUCTIONS

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL by Blue Force Gear. This medical kit is designed to provide essential supplies for emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the kit.

General Safety Guidelines

- **Product Safety:** This kit is designed for safe use in nonfood emergency situations. Ensure that all components are used as intended.
- **Enhanced Recalls:** In the event of a safety issue, updates and recall notices will be provided. Always check for the latest information.
- **Online Shopping:** If purchased online, this product is subject to the same safety regulations as items bought in physical stores.
- **Special Consumer Focus:** Extra care should be taken when using this kit around vulnerable groups, including children.
- **EU Contact Point:** For any safety inquiries, please refer to the product information provided with your purchase.
- **Rapid Alerts:** Stay informed about product safety updates through the EU Safety Gate system.

Specific Safety Precautions for Use

- **Age Restrictions:** This kit is not intended for use by children without adult supervision.
- **Check Contents:** Before use, check that all items are present and undamaged. Do not use any item that appears compromised.
- **Proper Use:** Use items in the kit only for their intended purpose. For example, do not use medical dressings for nonmedical purposes.
- **Hygiene:** Always use gloves when handling open wounds or bodily fluids to minimize the risk of infection.
- **Training:** It is recommended that users have basic first aid training to effectively utilize the contents of the kit.

Instructions for Installation and Usage

- **Accessing the Kit:**
 - Open the MOLLE pouch by pulling the top flap to access the contents.
 - Ensure you are in a safe environment when retrieving items.
- **Using the Contents:**
 - **Trauma Gloves:** Put on gloves before treating any injuries.
 - **QuikClot Bleeding Control Dressing:** Apply directly to the wound and apply pressure as needed.
 - **Emergency Trauma Dressing:** Use to cover larger wounds, ensuring it is secured in place.
 - **Hypothermia Blanket:** Wrap the injured person to maintain body temperature.
 - **Flat Fold Tape:** Use to secure dressings or gauze in place.
 - **Compressed Gauze:** Pack into wounds to control bleeding.
 - **Full Size Chest Seal:** Apply to chest wounds to prevent air from entering the chest cavity.
 - **Nasopharyngeal Airway:** Insert carefully to maintain an open airway if necessary.
 - **Decompression Needle:** Only use if trained; this is for advanced use in cases of tension

pneumothorax.

- **After Use:**

- Dispose of used items according to local regulations.
- Replenish any used supplies as soon as possible to maintain the kit's readiness.

Disposal Instructions

- Dispose of all used medical supplies in accordance with local health regulations.
- Ensure that all items are sealed in a bag before disposal to prevent contamination.
- Follow local guidelines for the disposal of hazardous materials, such as sharps.

Contact Information for Further Support

For any further inquiries or support regarding the TRAUMA KIT NOW! SMALL, please refer to the product documentation included with your purchase.

Thank you for your attention to these safety instructions. Proper use of the TRAUMA KIT NOW! SMALL can help ensure safety and effectiveness in emergency situations.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk