

TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL-MOLLE- ADVANCED SUPPLIES-WOLF GRAY

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL-MOLLE- ADVANCED SUPPLIES-WOLF GRAY](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101030
- Mfr. No.: HW-TKN-SMALL-AD
- Size: 1 Person
- Delivery weight: 0.442kg
- Shipping height: 102mm
- Shipping width: 178mm
- Shipping length: 254mm
- UPC: 810073653351

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [TRAUMA KIT NOW! SMALL BLUE FORCE GEAR TRAUMA KIT NOW! SMALLMOLLE ADVANCED SUPPLIESWOLF GRAY Safety Instruction Guide](#)
- [About Us](#)

TRAUMA KIT NOW! SMALL BLUE FORCE GEAR TRAUMA KIT NOW! SMALLMOLLE ADVANCED SUPPLIESWOLF GRAY Safety Instruction Guide

Introduction

The TRAUMA KIT NOW! SMALL is an advanced medical kit designed for emergency situations. It is manufactured by Blue Force Gear and is specifically tailored to provide essential medical supplies in a compact format. This guide outlines important safety instructions and guidelines to ensure safe and effective use of the kit, in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- Ensure that the trauma kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace any expired items immediately.
- Keep the kit out of reach of children unless under adult supervision.
- Always read and understand the instructions before using any medical supplies included in the kit.
- In case of an emergency, call for professional medical help as needed.

Specific Safety Precautions for Use

- Trauma Gloves:
 - Use gloves to prevent contamination when treating wounds.
 - Dispose of gloves properly after use to avoid crosscontamination.
- QuikClot Bleeding Control Dressing:
 - Apply pressure to the wound after placing the dressing to control bleeding effectively.
 - Do not apply to severe burns or deep puncture wounds without professional assistance.
- Emergency Trauma Dressing:
 - Ensure the dressing is applied tightly to control bleeding.
 - Seek professional medical help after applying the dressing.
- Hypothermia Blanket:
 - Use the blanket to maintain body temperature in cases of shock or exposure to cold.
 - Ensure the person is in a safe position before applying the blanket.
- Flat Fold Tape and Compressed Gauze:
 - Use tape to secure dressings in place.
 - Always use compressed gauze for bleeding control; do not use expired gauze.

Instructions for Installation and Usage

1. Initial Setup:

- Open the trauma kit and familiarize yourself with the contents.
- Ensure all items are present and in good condition.

2. Using the Kit:

- In an emergency, assess the situation quickly and ensure your safety first.
- Put on trauma gloves before handling any medical supplies.
- If there is bleeding:
 - Use QuikClot or compressed gauze to apply direct pressure to the wound.
 - If bleeding does not stop, apply additional dressings without removing the first one.
- For traumatic injuries:

- Use the emergency trauma dressing as needed.
- Secure it tightly but not so tight as to cut off circulation.
- If the person is at risk of hypothermia:
 - Wrap them in the hypothermia blanket to preserve body heat.
- After providing care, dispose of used materials safely and wash your hands thoroughly.

3. **PostUse Care:**

- Replenish any used items in the kit as soon as possible.
- Store the kit in its designated location for quick access in future emergencies.

Disposal Instructions

- Dispose of used medical supplies, including gloves and dressings, in accordance with local regulations for hazardous waste.
- Do not dispose of medical supplies in regular household waste.
- Ensure that any used items are sealed in a plastic bag before disposal to prevent contamination.

Contact Information for Further Support

For any inquiries or concerns regarding the TRAUMA KIT NOW! SMALL, please refer to the contact information provided with your product packaging or visit the manufacturer's website for additional resources and support.

By following these safety instructions and guidelines, you can ensure the effective use of the TRAUMA KIT NOW! SMALL while adhering to safety regulations. Always prioritize safety and seek professional medical assistance when necessary.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk