## **INSERT - MultiCam®**

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# INSERT MultiCam® Safety Instructions

## Introduction

Thank you for choosing the INSERT MultiCam® Individual First Aid Kit. This guide provides essential safety instructions to ensure the safe and effective use of your first aid kit. Please read this manual carefully before use.

## General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace any expired items.
- Keep the kit out of reach of children unless supervised by an adult.
- Use the kit only for its intended purpose and as described in this manual.
- In case of any doubt regarding the use of any item in the kit, consult a medical professional.

## Specific Safety Precautions for Use

- Always wash your hands before using any medical supplies from the kit to prevent infection.
- Avoid using any items that appear damaged or expired.
- Be aware of any allergies you or others may have to specific medical supplies (e.g., latex gloves).
- Follow the MARCH protocol (Massive hemorrhage, Airway, Respiratory, Circulation, Head injury) when providing first aid.
- If you are unsure about how to use any item in the kit, seek guidance from a trained professional.

## Instructions for Installation and Usage

### 1. Unpacking the Kit

- Carefully open the kit and check that all items listed in the contents are present.
- Inspect each item for any signs of damage or expiration.

### 2. Using the Kit

- In case of an emergency, quickly assess the situation and determine the necessary first aid steps.
- Follow the MARCH protocol:
  - **Massive Hemorrhage:** Apply direct pressure to control bleeding.
  - **Airway:** Ensure the airway is clear.
  - **Respiratory:** Check for breathing and provide assistance if necessary.
  - **Circulation:** Monitor pulse and circulation.
  - **Head Injury:** Be cautious with head injuries and seek professional help if needed.

### 3. After Use

- Dispose of any used items (e.g., bandages, gloves) in accordance with local regulations.
- Replenish any used supplies as soon as possible to maintain the kit's readiness.

## Disposal Instructions

- Dispose of medical waste according to local regulations. This may include:
  - Used bandages and dressings.
  - Expired medications.
  - Any other items that may be considered hazardous.
- Do not dispose of medical waste in regular household trash unless permitted by local laws.

## Contact Information for Further Support

For safety inquiries or further assistance regarding the INSERT MultiCam® kit, please refer to the EUbased contact point provided with your product.

Thank you for taking the time to read these safety instructions. Your safety and wellbeing are our top priority.

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