

INSERT - MultiCam®

INSERT è il nuovo kit individuale di primo soccorso progettato per contenere i presidi medici salvavita, necessari all'attuazione del protocollo MARCH. Sviluppato come inserto interno del kit RIFAK ([Horizontal](#) - [Vertical](#)), può essere acquistato separatamente per un utilizzo "stand-alone".



Table of Contents

- [Startpage](#)
- [INSERT MultiCam® Safety Instructions](#)
- [About Us](#)

INSERT MultiCam® Safety Instructions

Introduction

Thank you for choosing the INSERT MultiCam® Individual First Aid Kit. This guide provides essential safety instructions to ensure the safe and effective use of your first aid kit. Please read this manual carefully before use.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit for expiration dates and replace any expired items.
- Keep the kit out of reach of children unless supervised by an adult.
- Use the kit only for its intended purpose and as described in this manual.
- In case of any doubt regarding the use of any item in the kit, consult a medical professional.

Specific Safety Precautions for Use

- Always wash your hands before using any medical supplies from the kit to prevent infection.
- Avoid using any items that appear damaged or expired.
- Be aware of any allergies you or others may have to specific medical supplies (e.g., latex gloves).
- Follow the MARCH protocol (Massive hemorrhage, Airway, Respiratory, Circulation, Head injury) when providing first aid.
- If you are unsure about how to use any item in the kit, seek guidance from a trained professional.

Instructions for Installation and Usage

1. Unpacking the Kit

- Carefully open the kit and check that all items listed in the contents are present.
- Inspect each item for any signs of damage or expiration.

2. Using the Kit

- In case of an emergency, quickly assess the situation and determine the necessary first aid steps.
- Follow the MARCH protocol:
 - **Massive Hemorrhage:** Apply direct pressure to control bleeding.
 - **Airway:** Ensure the airway is clear.
 - **Respiratory:** Check for breathing and provide assistance if necessary.
 - **Circulation:** Monitor pulse and circulation.
 - **Head Injury:** Be cautious with head injuries and seek professional help if needed.

3. After Use

- Dispose of any used items (e.g., bandages, gloves) in accordance with local regulations.
- Replenish any used supplies as soon as possible to maintain the kit's readiness.

Disposal Instructions

- Dispose of medical waste according to local regulations. This may include:
 - Used bandages and dressings.
 - Expired medications.

- Any other items that may be considered hazardous.
- Do not dispose of medical waste in regular household trash unless permitted by local laws.

Contact Information for Further Support

For safety inquiries or further assistance regarding the INSERT MultiCam® kit, please refer to the EUbased contact point provided with your product.

Thank you for taking the time to read these safety instructions. Your safety and wellbeing are our top priority.

About Us

Brownells Italia

Brownells Italia - The shop of the largest dealer of components and accessories for target shooting guns, defense and hunting

**Sede Legale e Operativa **

Brownells Italia Srl

Via Dante Alighieri 163

25068 Sarezzo Bs

P.Iva e C.F. IT02801420981

www.brownells.it