

RIFAK Horizontal - Coyote

RIFAK (Rip-Away IFAK) è il nuovo kit individuale di primo soccorso progettato per contenere i presidi medici salvavita necessari all'attuazione del protocollo MARCH. La RIFAK conta due componenti distinte: un involucro esterno e un inserto interno estraibile, dove sono contenuti i presidi medici. Il sistema di sgancio è stato sviluppato per essere estremamente sicuro ma anche intuitivo da estrarre.



CARATTERISTICHE:

- Grazie al sistema Rip-Away l'operatore sarà in grado di raggiungere tutti i presidi medici in meno di 6 secondi.
- Non c'è rischio che i presidi medici possano fuoriuscire accidentalmente dall'inserto, anche estraendolo bruscamente, perché sono vincolati su tutti i lati.
- Non sono stati aggiunti "abbellimenti" (es. croci mediche, velcro, ricami...) che possano confondere l'operatore. L'unica porzione evidenziata in rosso è la maniglia di estrazione.
- Il cordino elastico esterno può essere usato per comprimere le dimensioni della tasca oltre che per trasportare un Tourniquet.
- Sui PALS/MOLLE laterali è possibile installare ulteriori tasche e/o vincolare strumenti come Forbici, Strap-Cutter e Tourniquet.

Table of Contents

- [Startpage](#)
- [RIFAK Horizontal Coyote Safety Instructions](#)
- [About Us](#)

RIFAK Horizontal Coyote Safety Instructions

Introduction

Thank you for choosing the RIFAK Horizontal Coyote, an innovative individual first aid kit designed to provide essential medical supplies for emergency situations. This guide outlines important safety instructions to ensure the proper use and maintenance of your RIFAK kit, in compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- **Product Safety:** Ensure safe use of the RIFAK kit by following all instructions provided in this manual.
- **Enhanced Recalls:** Stay informed about product recalls. If you receive a recall notice, follow the instructions provided to return or replace the product.
- **Online Shopping:** When purchasing online, ensure that the seller meets all safety requirements.
- **Special Consumer Focus:** This product is not intended for children. Keep out of reach of minors.
- **EU Contact Point:** For safety inquiries, contact the designated EU contact point associated with this product.
- **Rapid Alerts:** Monitor the EU's Safety Gate platform for updates on product safety.

Specific Safety Precautions for Use

- **Avoid Overloading:** Do not exceed the recommended capacity of the kit. Overloading may compromise the integrity of the components.
- **Inspect Regularly:** Check the kit regularly for any signs of damage or wear. Replace any items that are expired or damaged.
- **Proper Handling:** Handle the kit with care. Avoid dropping or exposing it to extreme temperatures or moisture.
- **Training Recommended:** It is highly recommended that users receive training in first aid and emergency response to effectively utilize the kit.
- **Use as Intended:** The RIFAK kit is designed for emergency first aid. Do not use it for nonmedical purposes.

Instructions for Installation and Usage

1. Unpacking the Kit:

- Carefully remove the RIFAK kit from its packaging.
- Inspect the contents to ensure all components are present and undamaged.

2. Accessing Medical Supplies:

- Use the RipAway system to quickly access medical supplies in under 6 seconds.
- Pull the red extraction handle to remove the inner insert containing the medical supplies.

3. Using the Medical Supplies:

- Follow the instructions provided with each medical supply for proper usage.
- Use supplies only for their intended purposes.

4. Reinserting the Insert:

- After use, securely place the insert back into the outer casing.
- Ensure it is properly locked in to prevent accidental removal.

5. **External Cord Usage:**

- The external elastic cord can be used to compress the size of the pocket or to transport a Tourniquet.
- Ensure the cord is securely fastened during transport.

6. **Attachment of Additional Pockets:**

- Utilize the PALS/MOLLE system on the sides of the kit to attach additional pockets or secure tools like scissors, strapcutters, and tourniquets.
- Ensure attachments are secure to prevent loss during use.

Disposal Instructions

- **Disposal of Expired Supplies:** Dispose of any expired or damaged medical supplies in accordance with local regulations.
- **Kit Disposal:** If the kit is no longer usable, dispose of it responsibly. Check with local waste management authorities for proper disposal methods.
- **Environmental Considerations:** Recycle materials where possible.

Contact Information for Further Support

For any inquiries regarding the RIFAK Horizontal Coyote, please refer to the contact information provided by the retailer or manufacturer. Ensure you have your purchase details ready for assistance.

By following these safety instructions, you can ensure the effective and safe use of your RIFAK Horizontal Coyote first aid kit. Stay prepared and safe.

About Us

Brownells Italia

Brownells Italia - The shop of the largest dealer of components and accessories for target shooting guns, defense and hunting

**Sede Legale e Operativa **

Brownells Italia Srl

Via Dante Alighieri 163

25068 Sarezzo Bs

P.Iva e C.F. IT02801420981

www.brownells.it