

# TRAUMA KIT NOW! ESSENTIALS SUPPLIES - BLUE FORCE GEAR

## TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNTED COYOTE

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - MEDIUM is a MOLLE mounted medical kit that can be deployed in seconds. This first aid pouch was originally designed to hold the contents of the Army IFAK and other medical supplies. This medical kit mounts on a 3 x 3 MOLLE field, a battle belt or other MOLLE Compatible Surfaces using the patented BFG Helium Whisper Attachment System. The Helium Whisper Attachment System is the lightest MOLLE attachment system in existence. Featuring less bulk so it sits closer to the body armor, pack or belt system. It has fewer parts to fail and is metal free to reduce secondary fragmentation. Best of all, Helium Whisper offers overall weight savings of up to 40% over traditional pouches. ESSENTIALS SUPPLIES INCLUDE: Trauma Gloves (Rolled 1x Pair) Quikclot Bleeding Control Dressing (1x) Emergency Trauma Dressing 6" (1x) Hypothermia Blanket - Green / Silver (1x) Flat Fold Tape (1x) Compressed Gauze (1x) SWAT-T Pressure Device (1x)



### Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! ESSENTIAL SUPPLIES MOLLE MOUNTED COYOTE](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 100050359
- Mfr. No.:
- Size: 1 Person
- Delivery weight: 0.635kg
- Shipping height: 76mm
- Shipping width: 178mm
- Shipping length: 229mm
- UPC: 810073650213

### Item details

Made in USA

# Table of Contents

- [Startpage](#)
- [Trauma Kit NOW! Essentials Supplies Safety Instruction Guide](#)
- [About Us](#)

# Trauma Kit NOW! Essentials Supplies Safety Instruction Guide

## Introduction

Thank you for purchasing the Trauma Kit NOW! – MEDIUM by Blue Force Gear. This MOLLE mounted medical kit is designed for rapid deployment and effective emergency response. This guide provides essential safety instructions to ensure safe use and compliance with the EU General Product Safety Regulation (GPSR).

## General Safety Guidelines

- **Product Safety:** Always ensure that the kit is used as intended. Familiarize yourself with the contents and their purposes.
- **Enhanced Recalls:** Stay informed about any product recalls. Check the EU Safety Gate platform regularly for updates.
- **Online Shopping:** If purchased online, ensure that the seller complies with safety regulations.
- **Special Consumer Focus:** This kit is not intended for use by children. Keep out of reach of minors.
- **EU Contact Point:** For safety inquiries, refer to the information provided by your retailer.
- **Rapid Alerts:** Be proactive in checking for alerts on unsafe products.

## Specific Safety Precautions for Use

- **Personal Safety:** Always wear protective gloves when handling medical supplies to prevent contamination.
- **Storage:** Store the kit in a cool, dry place, away from direct sunlight and moisture to maintain the integrity of the supplies.
- **Inspection:** Regularly check the kit for expired items or damage. Replace any used or expired components immediately.
- **Training:** Ensure that users are trained in basic first aid and CPR techniques before using the kit in an emergency.

## Instructions for Installation and Usage

### Installation

- **MOLLE Attachment:**
  - Locate a 3 x 3 MOLLE field on your tactical gear, battle belt, or compatible surface.
  - Align the Trauma Kit NOW! with the MOLLE webbing.
  - Use the patented BFG Helium Whisper Attachment System to secure the kit. Ensure it is firmly attached to prevent accidental detachment.

### Usage

- **Accessing the Kit:**
  - In an emergency, quickly release the kit from its mounted position.
  - Open the pouch to access the supplies.
- **Contents:**
  - **Trauma Gloves:** Use these gloves to protect yourself and the patient.
  - **Quikclot Bleeding Control Dressing:** Apply to severe bleeding areas.
  - **Emergency Trauma Dressing:** Use for major wounds.

- **Hypothermia Blanket:** Keep the patient warm if necessary.
- **Flat Fold Tape:** Secure dressings in place.
- **Compressed Gauze:** Use for wound packing.
- **SWATT Pressure Device:** Employ for effective bleeding control.

- **PostUse:**

- Dispose of used items safely and in accordance with local regulations.
- Restock the kit after use to ensure it is ready for future emergencies.

## **Disposal Instructions**

- Dispose of all used medical supplies in accordance with local regulations.
- Follow guidelines for the disposal of hazardous materials if applicable.
- Ensure that the kit itself is kept intact for future use.

## **Contact Information for Further Support**

For any inquiries or support regarding the Trauma Kit NOW! – MEDIUM, please refer to the contact information provided by your retailer. Regularly check for updates on safety and product recalls through the EU Safety Gate platform.

By following these guidelines, you can ensure safe and effective use of your Trauma Kit NOW! – MEDIUM, while adhering to EU safety regulations. Thank you for prioritizing safety in your emergency preparedness efforts.

# About Us

## Brownells UK

### **Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories**

Unit 1, Laughing Dog Industrial Estate  
London Road  
Rugby  
Warwickshire  
CV23 9LP

[www.brownells.co.uk](http://www.brownells.co.uk)