

TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - PRO SUPPLIES-MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - PRO SUPPLIES-MULTICAM](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101023
- Mfr. No.: HW-TKN-SMALL-PRO-MC
- Size: 1 Person
- Delivery weight: 0.44kg
- Shipping height: 51mm
- Shipping width: 203mm
- Shipping length: 279mm
- UPC: 810073653283

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Safety Instruction Guide for TRAUMA KIT NOW! SMALL](#)
- [About Us](#)

Safety Instruction Guide for TRAUMA KIT NOW! SMALL

Introduction

Thank you for choosing the TRAUMA KIT NOW! SMALL from Blue Force Gear. This medical pouch is designed to provide essential first aid supplies in a compact format. It is important to understand how to use this kit safely and effectively. This guide provides essential safety instructions to ensure your safety and the safety of others.

General Safety Guidelines

- **Product Safety:** Ensure safe use of the Trauma Kit NOW! SMALL by following the instructions provided in this guide. Always check the contents for any damage before use.
- **Enhanced Recalls:** Stay informed about product recalls. If you receive a recall notice, follow the instructions provided to return or exchange the product.
- **Online Shopping:** When purchasing online, ensure that the seller complies with safety requirements and provides clear product information.
- **Special Consumer Focus:** This kit is not intended for children. Keep out of reach of minors to prevent accidental use.
- **EU Contact Point:** For any safety inquiries, please refer to your local EU consumer protection agency.
- **Rapid Alerts:** Stay updated on unsafe products through the EU Safety Gate platform. Regularly check for any alerts regarding the Trauma Kit NOW! SMALL.

Specific Safety Precautions for Use

- **Read the Instructions:** Before using the kit, read all instructions carefully to understand how to use the contents safely.
- **Inspect Contents:** Regularly check the kit for expired or damaged items. Replace any items as necessary to ensure the kit is ready for use.
- **Hygiene Practices:** Always use clean hands when handling medical supplies. If possible, wear gloves to prevent contamination.
- **Emergency Situations:** In case of a medical emergency, use the kit as directed. If the situation is lifethreatening, seek professional medical assistance immediately.
- **Storage:** Store the kit in a cool, dry place away from direct sunlight and extreme temperatures. Ensure it is easily accessible in case of an emergency.

Instructions for Installation and Usage

- **Installation:** The Trauma Kit NOW! SMALL can be attached to MOLLEcompatible gear. To install:
 - Locate the MOLLE webbing on your gear.
 - Slide the straps of the kit through the webbing.
 - Secure the kit by pulling the straps tight to ensure it does not move during use.
- **Usage:** To access the contents of the kit:
 - Pull the top of the pouch in a single motion to open.
 - Remove the necessary items carefully to avoid spilling the contents.
 - Follow the instructions provided with each item for proper use.

Disposal Instructions

- **Expired Items:** Dispose of expired or damaged items in accordance with local regulations. Do not throw them in regular trash if they contain hazardous materials.
- **Kit Disposal:** If the entire kit is no longer needed, check with local waste management services for proper disposal methods. Ensure that any medical supplies are disposed of safely to prevent harm to

others.

Contact Information for Further Support

For further assistance or inquiries regarding the TRAUMA KIT NOW! SMALL, please refer to your local consumer protection agency or the retailer from whom you purchased the kit.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk