

TRAUMA KIT NOW! - SMALL - BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - ADVANCED SUPPLIES-MULTICAM

[Emergency & Survival Gear](#) > [First Aid](#) > [First Aid Kits](#)

The Trauma Kit NOW! - SMALL is the newest med kit from Blue Force Gear. Designed in response to customers who loved the features and design of the TKN-MEDIUM but needed a smaller pouch. This MOLLE medical pouch is 30% smaller than the TKN-MEDIUM. The TKN-SMALL utilizes the single-motion pull out the top design of the original TKN but in the familiar size and format of a Double / Triple (stacked) M4 magazine pouch. ESSENTIALS INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Bleeding Control Dressing (1x - white), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x) Flat Fold Tape (1x), and Compressed Gauze (1x) PR INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze LE (1x - black), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), and Full Size Chest Seal (1x Twin Pack) ADVANCED INCLUDES: Trauma Gloves (Rolled 2x Pair), QuikClot Combat Gauze (1x - green), Emergency Trauma Dressing 6" (1x), Hypothermia Blanket (1x), Flat Fold Tape (1x), Compressed Gauze (1x), Full Size Chest Seal (1x Twin Pack), Nasopharyngeal Airway (1x), and Decompression Needle (1x)



Attributes

- Name: [BLUE FORCE GEAR TRAUMA KIT NOW! SMALL - MOLLE - ADVANCED SUPPLIES-MULTICAM](#)
- Manufacturer: [BLUE FORCE GEAR](#)
- Product no.: 430101028
- Mfr. No.: HW-TKN-SMALL-ADV-MC
- Size: 1 Person
- Delivery weight: 0.567kg
- Shipping height: 279mm
- Shipping width: 76mm
- Shipping length: 229mm
- UPC: 810073653337

Item details

Made in USA

Table of Contents

- Startpage
- Safety Instruction Guide for the Trauma Kit NOW! SMALL
- About Us

Safety Instruction Guide for the Trauma Kit NOW! SMALL

Introduction

Thank you for choosing the Trauma Kit NOW! SMALL from Blue Force Gear. This compact medical pouch is designed to provide essential trauma supplies in a convenient size. It is important to understand how to use this product safely and effectively. This guide outlines safety guidelines, precautions, and instructions to ensure your safety and the safety of others.

General Safety Guidelines

- Always read and follow the instructions provided in this guide before using the Trauma Kit NOW! SMALL.
- Ensure that the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure that all items are present and within their expiration dates.
- Keep the kit out of reach of children unless under adult supervision.
- In case of an emergency, always seek professional medical assistance in addition to using the kit.
- Report any unsafe products or accidents to the relevant authorities.
- Stay informed about product recalls by checking the EU's Safety Gate platform.

Specific Safety Precautions for Use

- **Trauma Gloves:** Always wear the trauma gloves when handling blood or bodily fluids to prevent contamination and infection.
- **QuikClot Products:** Follow the instructions on the packaging carefully to ensure proper application for bleeding control.
- **Emergency Trauma Dressing:** Use only as directed. If bleeding does not stop, seek immediate medical assistance.
- **Hypothermia Blanket:** Ensure the blanket is used correctly to provide warmth in emergency situations.
- **Compressed Gauze:** Use as directed for wound packing and control of bleeding.
- **Chest Seal:** Ensure proper application to prevent air from entering the chest cavity in case of a penetrating chest injury.
- **Nasopharyngeal Airway and Decompression Needle:** These items should only be used by trained personnel. Improper use can lead to serious injury.

Instructions for Installation and Usage

1. Opening the Kit:

- The Trauma Kit NOW! SMALL features a single motion pullout top design.
- To access the contents, pull the top flap upwards smoothly.

2. Using the Supplies:

- **Trauma Gloves:** Roll out the gloves from the pack and wear them before treating any wounds.
- **QuikClot Bleeding Control Dressing:** Open the packaging and apply directly to the wound, applying pressure as necessary.
- **Emergency Trauma Dressing:** Remove from packaging and apply over the wound, securing it in place.
- **Hypothermia Blanket:** Unfold and wrap around the individual needing warmth.
- **Compressed Gauze:** Pack tightly into the wound to control bleeding.
- **Chest Seal:** Remove the protective backing and apply over the wound, ensuring a tight seal.

- **Nasopharyngeal Airway:** Lubricate and gently insert into the nostril until it reaches the back of the throat.
- **Decompression Needle:** Only trained personnel should perform this procedure. Locate the site and insert the needle to relieve pressure.

3. PostUse:

- Dispose of used gloves and medical supplies in accordance with local regulations.
- Clean the kit and store it properly for future use.

Disposal Instructions

- Dispose of any expired or used items in accordance with local waste disposal regulations.
- Do not dispose of medical supplies in regular trash unless specifically instructed to do so.
- Follow local guidelines for the disposal of hazardous materials when applicable.

Contact Information for Further Support

For any questions, concerns, or further assistance regarding the Trauma Kit NOW! SMALL, please refer to the contact information provided with your product packaging or the manufacturer's website.

Thank you for prioritizing safety and preparedness with the Trauma Kit NOW! SMALL.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk