

DELTA BAG - SPIRITUS SYSTEMS DELTA BAG WOODLAND

[Emergency & Survival Gear](#) > [Backpacks & Bags](#) > [Backpacks](#)

The Delta Bag is a compact, assault medic bag that is designed to be worn over a plate carrier and hold a wide variety of medical supplies and equipment. The front of the bag is divided into three uniformly sized compartments with diagonal zipper access that enables the end user to quickly access the entire pocket without having to reposition any of its contents. The unique configuration also aids in passive retention of gear when the compartment is opened. The main compartment of the Delta Bag is a clamshell design with two full fields of loop that can accept any hookbacked accessories (including the entire line of Spiritus Systems utility inserts, and other 3rd party products). The bag can be transported with the clamshell in the open position from location to location using the integrated briefcase handles, and provides an intuitive workstation when emergency medical care is being administered. The back side of the pack features two slim compartments for keeping flat-pack items secure and damage free, and keeps the pack close to the body when worn. The included straps for the delta bag are removeable and replaceable with any 3rd party straps that attach with 1" webbing. The bag features four sewn-in metal tri-glide mounting points for strap integration. The Delta bag features mesh and metal eyelets in key locations for full water drainage capability, and is designed to remain effective and accessible in austere environments.



Attributes

- Name: [SPIRITUS SYSTEMS DELTA BAG WOODLAND](#)
- Manufacturer: [SPIRITUS SYSTEMS](#)
- Product no.: 430104969
- Mfr. No.: BU-DB-WD
- Colour: Woodland
- Weight: 1.261 lbs
- Delivery weight: 0.587kg
- Shipping height: 76mm
- Shipping width: 457mm
- Shipping length: 305mm

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [Delta Bag Safety Instruction Guide](#)
- [About Us](#)

Delta Bag Safety Instruction Guide

Introduction

Thank you for choosing the Delta Bag by Spiritus Systems. This compact assault medic bag is designed for efficient and effective use in emergency medical situations. To ensure your safety and the longevity of the product, please read this safety instruction guide carefully.

General Safety Guidelines

- Ensure that the Delta Bag is used only for its intended purpose as a medic bag.
- Regularly inspect the bag for any signs of wear or damage before each use.
- Keep the bag out of reach of children when not in use.
- Always store the bag in a cool, dry place to prevent damage from moisture or extreme temperatures.
- Follow all local regulations regarding the transport and use of medical supplies.

Specific Safety Precautions for Use

- When accessing the compartments, ensure the bag is securely positioned to prevent it from tipping over.
- Be cautious when handling sharp objects or medical equipment stored in the bag.
- Avoid overloading the bag beyond its intended capacity, as this may compromise its integrity and safety.
- When using the bag in austere environments, be aware of your surroundings to avoid potential hazards.
- If the bag becomes wet, allow it to dry completely before storing to prevent mold and mildew.

Instructions for Installation and Usage

• Attaching Straps:

- Remove the existing straps if necessary.
- Align the 1" webbing of the new strap with the metal triglide mounting points on the bag.
- Securely fasten the strap to ensure it can support the weight of the bag when in use.

• Accessing Compartments:

- Open the front diagonal zippers to access the three compartments.
- Ensure that items within the compartments are organized and secured to prevent movement during transport.

• Using the Clamshell Design:

- To open the main compartment, unzip the clamshell design fully.
- Utilize the loop fields to attach any hookbacked accessories as needed.
- When transporting, you may carry the bag with the clamshell open using the integrated briefcase handles.

• Maintaining the Bag:

- Regularly clean the exterior and interior of the bag with a damp cloth and mild detergent.
- Inspect the mesh and metal eyelets for any blockages to ensure proper drainage.

Disposal Instructions

- When the Delta Bag reaches the end of its useful life, dispose of it responsibly.
- If the bag is damaged beyond repair, consider recycling the materials where possible.

- Follow local regulations for the disposal of medical supplies and equipment.

Contact Information for Further Support

For any safety inquiries, product support, or to report unsafe products, please refer to the official EU Safety Gate platform or your local consumer safety authority.

About Us

Brownells UK

Brownells UK - World's Largest Supplier of Gun Parts, Gunsmith Tools & Shooting Accessories

Unit 1, Laughing Dog Industrial Estate
London Road
Rugby
Warwickshire
CV23 9LP

www.brownells.co.uk