

RIFAK Vertical - Black

RIFAK (Rip-Away IFAK) è il nuovo kit individuale di primo soccorso progettato per contenere i presidi medici salvavita necessari all'attuazione del protocollo MARCH. La RIFAK conta due componenti distinte: un involucro esterno e un inserto interno estraibile, dove sono contenuti i presidi medici. Il sistema di sgancio è stato sviluppato per essere estremamente sicuro ma anche intuitivo da estrarre.

CARATTERISTICHE:

- Grazie al sistema Rip-Away l'operatore sarà in grado di raggiungere tutti i presidi medici in meno di 6 secondi.
- Non c'è rischio che i presidi medici possano fuoriuscire accidentalmente dall'inserto, anche estraendolo bruscamente, perché sono vincolati su tutti i lati.
- Non sono stati aggiunti "abbellimenti" (es. croci mediche, velcro, ricami...) che possano confondere l'operatore. L'unica porzione evidenziata in rosso è la maniglia di estrazione.
- Il cordino elastico esterno può essere usato per comprimere le dimensioni della tasca oltre che per trasportare un Tourniquet.
- Sui PALS/MOLLE laterali è possibile installare ulteriori tasche e/o vincolare strumenti come Forbici, Strap-Cutter e Tourniquet.

Table of Contents

- [Startpage](#)
- [English: RIFAK Vertical Black Safety Instructions](#)
- [Italiano: Guida alle Istruzioni di Sicurezza del Prodotto](#)

RIFAK Vertical Black Safety Instructions

Introduction

Thank you for choosing the RIFAK Vertical Black, the innovative RipAway Individual First Aid Kit (IFAK) designed to provide essential medical supplies for emergency situations. This safety instruction guide is designed to help you understand how to use the product safely and effectively, ensuring compliance with the EU General Product Safety Regulation (GPSR).

General Safety Guidelines

- **Product Safety:** The RIFAK Vertical has been designed to ensure safe use. Please follow all instructions to minimize risks.
- **Enhanced Recalls:** In the unlikely event of a product recall, standardized notices will be provided. Please stay informed about recalls.
- **Online Shopping:** If purchased online, your product is subject to the same safety protections as those bought in physical stores.
- **Special Consumer Focus:** Extra care should be taken when using this product around vulnerable groups, including children.
- **EU Contact Point:** For any safety inquiries, please refer to the designated EU contact point.
- **Rapid Alerts:** Stay updated on product safety through the EU's Safety Gate system.

Specific Safety Precautions for Use

- Ensure the RIFAK Vertical is stored in a secure location to prevent unauthorized access, especially by children.
- Familiarize yourself with the contents and features of the kit before an emergency arises.
- Practice using the RipAway system to ensure you can access medical supplies quickly and efficiently.
- Inspect the kit regularly for any signs of damage or wear. Replace any compromised components immediately.
- Do not overload the kit beyond its intended capacity, as this may compromise its functionality and safety.

Instructions for Installation and Usage

Installation

1. **Mounting Options:** The RIFAK Vertical can be mounted using PALS/MOLLE webbing. Ensure it is securely fastened to prevent accidental detachment.
2. **Attachment:** Use the external elastic cord to compress the pocket size or to carry a Tourniquet as needed.
3. **Additional Storage:** You may install additional pockets or secure tools such as scissors or strapcutters on the PALS/MOLLE sides.

Usage

1. **Accessing Supplies:**
 - To access the medical supplies, pull the red extraction handle firmly.
 - The RipAway system allows for quick access in less than 6 seconds.
 - Ensure that all medical supplies are securely held within the insert to prevent accidental spillage.
2. **Handling Medical Supplies:**
 - Use all medical supplies according to their specific instructions.
 - Always check the expiration dates on medical supplies and replace them as necessary.
3. **PostUse:**
 - After use, restock the kit with fresh supplies immediately.
 - Clean any reusable components according to the manufacturer's guidelines.

Disposal Instructions

- Dispose of any expired or used medical supplies according to local regulations.
- The RIFAK Vertical itself can be disposed of in regular waste unless otherwise specified by local disposal guidelines.
- Consider recycling components where applicable to minimize environmental impact.

Contact Information for Further Support

For any questions or concerns regarding the RIFAK Vertical Black, please refer to the designated EU contact point for safety inquiries.

By following these guidelines, you can ensure the safe and effective use of your RIFAK Vertical Black, contributing to your preparedness in emergency situations. Thank you for prioritizing safety!

Guida alle Istruzioni di Sicurezza del Prodotto

Introduzione

Benvenuti nella guida alle istruzioni di sicurezza per il RIFAK Vertical Black. Questo kit individuale di pronto soccorso è progettato per fornire un accesso rapido e sicuro ai dispositivi medici salvavita. Si prega di leggere attentamente questa guida per garantire un uso sicuro e corretto del prodotto.

Linee Guida Generali di Sicurezza

- Assicurarsi che il prodotto sia utilizzato solo per lo scopo previsto.
- Tenere il kit lontano dalla portata dei bambini se non supervisionati.
- Controllare regolarmente il contenuto del kit e sostituire eventuali dispositivi scaduti o danneggiati.
- In caso di emergenza, seguire sempre le linee guida di pronto soccorso e contattare i servizi di emergenza se necessario.
- Segnalare eventuali prodotti difettosi o incidenti alle autorità competenti.

Precauzioni di Sicurezza Specifiche per l'Uso

- Prima di utilizzare il RIFAK, assicurarsi di essere formati nell'uso dei dispositivi medici inclusi.
- Non forzare l'inserito interno durante l'estrazione; seguire le istruzioni di sgancio.
- Evitare di esporre il kit a temperature estreme o umidità eccessiva.
- Non utilizzare il kit se danneggiato o se l'involucro esterno è compromesso.
- Non utilizzare il kit in situazioni non di emergenza per garantire la disponibilità dei dispositivi quando necessario.

Istruzioni per l'Installazione e l'Uso

1. Installazione del Kit:

- Posizionare il RIFAK in un luogo facilmente accessibile, preferibilmente in un'area di pronto soccorso designata.
- Assicurarsi che il kit sia fissato correttamente se montato su un supporto PALS/MOLLE.

2. Uso del Kit:

- Per estrarre l'inserito interno, afferrare la maniglia evidenziata in rosso.
- Utilizzare il sistema RipAway per accedere rapidamente ai dispositivi medici.
- Verificare che i dispositivi siano saldamente bloccati prima di procedere all'uso.

3. Controllo dei Contenuti:

- Controllare il contenuto del kit ogni sei mesi.
- Sostituire i dispositivi scaduti o danneggiati immediatamente.

Istruzioni per lo Smaltimento

- Smaltire i dispositivi medici scaduti o danneggiati secondo le normative locali per i rifiuti sanitari.
- Non gettare il kit o i suoi componenti nell'ambiente. Seguire le linee guida locali per il riciclaggio e lo smaltimento.

Informazioni di Contatto per Ulteriori Supporto

Per ulteriori informazioni o domande riguardanti la sicurezza del prodotto, si prega di contattare il produttore o consultare il sito web ufficiale. Assicurarsi di avere il numero di modello e altre informazioni pertinenti a portata di mano.

Si prega di seguire queste linee guida per garantire un uso sicuro e efficace del RIFAK Vertical Black. La sicurezza è la nostra priorità.