

## RIFAK Vertical - Black

RIFAK ( Rip-Away IFAK ) è il nuovo kit individuale di primo soccorso progettato per contenere i presidi medici salvavita necessari all'attuazione del protocollo MARCH. La RIFAK conta due componenti distinte: un involucro esterno e un inserto interno estraibile, dove sono contenuti i presidi medici. Il sistema di sgancio è stato sviluppato per essere estremamente sicuro ma anche intuitivo da estrarre.



### CARATTERISTICHE:

- Grazie al sistema Rip-Away l'operatore sarà in grado di raggiungere tutti i presidi medici in meno di 6 secondi.
- Non c'è rischio che i presidi medici possano fuoriuscire accidentalmente dall'inserto, anche estraendolo bruscamente, perché sono vincolati su tutti i lati.
- Non sono stati aggiunti "abbellimenti" (es. croci mediche, velcro, ricami...) che possano confondere l'operatore. L'unica porzione evidenziata in rosso è la maniglia di estrazione.
- Il cordino elastico esterno può essere usato per comprimere le dimensioni della tasca oltre che per trasportare un Tourniquet.
- Sui PALS/MOLLE laterali è possibile installare ulteriori tasche e/o vincolare strumenti come Forbici, Strap-Cutter e Tourniquet.

# Table of Contents

- [Startpage](#)
- [RIFAK Vertical Black Safety Instructions](#)
- [About Us](#)

# RIFAK Vertical Black Safety Instructions

## Introduction

Thank you for choosing the RIFAK Vertical Black, the innovative RipAway Individual First Aid Kit (IFAK) designed to provide essential medical supplies for emergency situations. This safety instruction guide is designed to help you understand how to use the product safely and effectively, ensuring compliance with the EU General Product Safety Regulation (GPSR).

## General Safety Guidelines

- **Product Safety:** The RIFAK Vertical has been designed to ensure safe use. Please follow all instructions to minimize risks.
- **Enhanced Recalls:** In the unlikely event of a product recall, standardized notices will be provided. Please stay informed about recalls.
- **Online Shopping:** If purchased online, your product is subject to the same safety protections as those bought in physical stores.
- **Special Consumer Focus:** Extra care should be taken when using this product around vulnerable groups, including children.
- **EU Contact Point:** For any safety inquiries, please refer to the designated EU contact point.
- **Rapid Alerts:** Stay updated on product safety through the EU's Safety Gate system.

## Specific Safety Precautions for Use

- Ensure the RIFAK Vertical is stored in a secure location to prevent unauthorized access, especially by children.
- Familiarize yourself with the contents and features of the kit before an emergency arises.
- Practice using the RipAway system to ensure you can access medical supplies quickly and efficiently.
- Inspect the kit regularly for any signs of damage or wear. Replace any compromised components immediately.
- Do not overload the kit beyond its intended capacity, as this may compromise its functionality and safety.

## Instructions for Installation and Usage

### Installation

1. **Mounting Options:** The RIFAK Vertical can be mounted using PALS/MOLLE webbing. Ensure it is securely fastened to prevent accidental detachment.
2. **Attachment:** Use the external elastic cord to compress the pocket size or to carry a Tourniquet as needed.
3. **Additional Storage:** You may install additional pockets or secure tools such as scissors or strapcutters on the PALS/MOLLE sides.

### Usage

1. **Accessing Supplies:**
  - To access the medical supplies, pull the red extraction handle firmly.
  - The RipAway system allows for quick access in less than 6 seconds.
  - Ensure that all medical supplies are securely held within the insert to prevent accidental spillage.
2. **Handling Medical Supplies:**
  - Use all medical supplies according to their specific instructions.

- Always check the expiration dates on medical supplies and replace them as necessary.

### 3. **PostUse:**

- After use, restock the kit with fresh supplies immediately.
- Clean any reusable components according to the manufacturer's guidelines.

## **Disposal Instructions**

- Dispose of any expired or used medical supplies according to local regulations.
- The RIFAK Vertical itself can be disposed of in regular waste unless otherwise specified by local disposal guidelines.
- Consider recycling components where applicable to minimize environmental impact.

## **Contact Information for Further Support**

For any questions or concerns regarding the RIFAK Vertical Black, please refer to the designated EU contact point for safety inquiries.

By following these guidelines, you can ensure the safe and effective use of your RIFAK Vertical Black, contributing to your preparedness in emergency situations. Thank you for prioritizing safety!

# About Us

## Brownells Italia

**Brownells Italia - The shop of the largest dealer of components and accessories for target shooting guns, defense and hunting**

**Sede Legale e&nbsp; Operativa&nbsp;**

Brownells Italia Srl&nbsp;

Via Dante Alighieri 163

25068&nbsp; Sarezzo Bs

P.Iva e C.F. IT02801420981

[www.brownells.it](http://www.brownells.it)