

RIFAK Vertical - Coyote

RIFAK (Rip-Away IFAK) è il nuovo kit individuale di primo soccorso progettato per contenere i presidi medici salvavita necessari all'attuazione del protocollo MARCH. La RIFAK conta due componenti distinte: un involucro esterno e un inserto interno estraibile, dove sono contenuti i presidi medici. Il sistema di sgancio è stato sviluppato per essere estremamente sicuro ma anche intuitivo da estrarre.



CARATTERISTICHE:

- Grazie al sistema Rip-Away l'operatore sarà in grado di raggiungere tutti i presidi medici in meno di 6 secondi.
- Non c'è rischio che i presidi medici possano fuoriuscire accidentalmente dall'inserto, anche estraendolo bruscamente, perché sono vincolati su tutti i lati.
- Non sono stati aggiunti "abbellimenti" (es. croci mediche, velcro, ricami...) che possano confondere l'operatore. L'unica porzione evidenziata in rosso è la maniglia di estrazione.
- Il cordino elastico esterno può essere usato per comprimere le dimensioni della tasca oltre che per trasportare un Tourniquet.
- Sui PALS/MOLLE laterali è possibile installare ulteriori tasche e/o vincolare strumenti come Forbici, Strap-Cutter e Tourniquet.

Table of Contents

- [Startpage](#)
- [English: RIFAK Vertical Coyote Safety Instruction Guide](#)
- [Italiano: RIFAK Vertical Coyote: Guida alle Istruzioni di Sicurezza](#)

RIFAK Vertical Coyote Safety Instruction Guide

Introduction

Thank you for choosing the RIFAK Vertical Coyote, your reliable individual first aid kit designed to provide essential medical supplies for emergency situations. This guide aims to ensure your safe use of the product, in compliance with the EU General Product Safety Regulation (GPSR). Please read this guide carefully before using the kit.

General Safety Guidelines

- **Product Safety:** Ensure safe use of the RIFAK Vertical by following all instructions and guidelines provided in this manual.
- **Enhanced Recalls:** Stay informed about any recalls or safety notices related to this product. Regularly check the EU Safety Gate platform for updates.
- **Online Shopping:** If purchased online, ensure that the seller complies with safety requirements as mandated for physical stores.
- **Special Consumer Focus:** This product is not intended for use by children. Keep it out of reach to avoid accidental use.
- **EU Contact Point:** For any safety inquiries, please refer to the contact details provided by the manufacturer.
- **Rapid Alerts:** Be aware that updates on unsafe products will be communicated through the Safety Gate system.

Specific Safety Precautions for Use

- **Avoid Exposure to Extreme Conditions:** Do not expose the RIFAK Vertical to extreme temperatures, moisture, or direct sunlight for prolonged periods.
- **Inspect Before Use:** Always check the kit for any signs of damage or missing components before use.
- **Secure Handling:** Ensure that the kit is securely fastened during transport to prevent accidental opening.
- **Proper Use of Medical Supplies:** Use the medical supplies only for their intended purpose and according to the instructions provided with each item.

Instructions for Installation and Usage

- **Accessing the Kit:**
 - Utilize the RipAway system to quickly access medical supplies in less than 6 seconds.
 - Ensure that the extraction handle, highlighted in red, is easily visible and accessible for rapid deployment.
- **Using the External Elastic Cord:**
 - The external elastic cord can be used to compress the size of the pocket and securely carry a Tourniquet.
- **Attaching Additional Items:**
 - Use the PALS/MOLLE lateral attachment points to install extra pouches or secure tools such as scissors, strapcutters, and tourniquets as needed.
- **Storing the Kit:**

- Store the RIFAK Vertical in a cool, dry place when not in use to maintain the integrity of the medical supplies.

- **Regular Maintenance:**

- Periodically check and restock the medical supplies to ensure they are complete and within their expiration dates.

Disposal Instructions

- **Disposing of the Kit:** When the RIFAK Vertical is no longer usable, dispose of it in accordance with local regulations regarding the disposal of medical supplies and materials.
- **Recycling:** If applicable, separate recyclable components and dispose of them in the appropriate recycling bins.

Contact Information for Further Support

For any questions or concerns regarding the RIFAK Vertical Coyote, please refer to the manufacturer's contact information provided with your product.

Thank you for your attention to these safety instructions. Your safety and wellbeing are paramount, and we hope that the RIFAK Vertical Coyote serves you well in emergency situations.

RIFAK Vertical Coyote: Guida alle Istruzioni di Sicurezza

Introduzione

Grazie per aver scelto il kit RIFAK Vertical Coyote. Questo kit è progettato per fornire un pronto soccorso efficace in situazioni di emergenza. È importante seguire attentamente le istruzioni di sicurezza per garantire un utilizzo sicuro e corretto del prodotto.

Linee Guida Generali di Sicurezza

- Assicurati di utilizzare il kit RIFAK solo per scopi previsti e in conformità con le istruzioni fornite.
- Controlla regolarmente il contenuto del kit per garantire che sia completo e in buone condizioni.
- Non utilizzare prodotti scaduti o danneggiati.
- Conserva il kit in un luogo facilmente accessibile e visibile.
- In caso di emergenza, contatta i servizi di emergenza locali oltre all'uso del kit.

Precauzioni di Sicurezza Specifiche per l'Uso

- Prima di utilizzare il kit, familiarizzati con il contenuto e le funzioni di ciascun presidio medico.
- Assicurati di avere una buona presa sull'involucro esterno durante l'estrazione dell'insero interno.
- Non forzare l'insero se non si estrae facilmente; controlla che non ci siano ostruzioni.
- Non utilizzare il kit in condizioni di scarsa visibilità o in situazioni di pericolo imminente.
- Proteggi il kit dall'umidità e dalle temperature estreme per preservare l'integrità dei materiali.

Istruzioni per l'Installazione e l'Uso

1. Installazione del Kit:

- Posiziona il kit in un luogo facilmente accessibile, come un veicolo o un'area di lavoro.
- Utilizza il cordino elastico esterno per fissare il kit a un supporto sicuro, se necessario.

2. Uso del Kit:

- In caso di emergenza, attiva il sistema RipAway per accedere rapidamente ai presidi medici.
- Estrarre l'insero interno con cautela, assicurandosi che tutti i materiali siano ben fissati.
- Utilizza i presidi medici secondo le istruzioni specifiche fornite con ciascun articolo.
- Se necessario, utilizza le tasche PALS/MOLLE per aggiungere strumenti o materiali supplementari.

Istruzioni per lo Smaltimento

- Smaltisci i materiali del kit secondo le normative locali riguardanti i rifiuti pericolosi.
- Non gettare i prodotti scaduti o danneggiati nei rifiuti domestici.
- Contatta un centro di smaltimento autorizzato per informazioni su come smaltire in modo sicuro i materiali medici.

Informazioni di Contatto per Ulteriori Supporto

Per domande o ulteriori informazioni sulla sicurezza del prodotto, ti invitiamo a contattare il tuo rivenditore o produttore. Assicurati di conservare la prova d'acquisto e di fare riferimento al numero di modello del prodotto per un'assistenza più rapida.

Si prega di notare che il kit RIFAK Vertical Coyote è progettato per l'uso da parte di operatori formati. La

sicurezza e l'efficacia del kit dipendono dall'uso corretto e dalla comprensione dei materiali in esso contenuti.