

RIFAK Vertical - Ranger Green

RIFAK (Rip-Away IFAK) è il nuovo kit individuale di primo soccorso progettato per contenere i presidi medici salvavita necessari all'attuazione del protocollo MARCH. La RIFAK conta due componenti distinte: un involucro esterno e un inserto interno estraibile, dove sono contenuti i presidi medici. Il sistema di sgancio è stato sviluppato per essere estremamente sicuro ma anche intuitivo da estrarre.



CARATTERISTICHE:

- Grazie al sistema Rip-Away l'operatore sarà in grado di raggiungere tutti i presidi medici in meno di 6 secondi.
- Non c'è rischio che i presidi medici possano fuoriuscire accidentalmente dall'inserto, anche estraendolo bruscamente, perché sono vincolati su tutti i lati.
- Non sono stati aggiunti "abbellimenti" (es. croci mediche, velcro, ricami...) che possano confondere l'operatore. L'unica porzione evidenziata in rosso è la maniglia di estrazione.
- Il cordino elastico esterno può essere usato per comprimere le dimensioni della tasca oltre che per trasportare un Tourniquet.
- Sui PALS/MOLLE laterali è possibile installare ulteriori tasche e/o vincolare strumenti come Forbici, Strap-Cutter e Tourniquet.

Table of Contents

- [Startpage](#)
- [RIFAK Vertical Ranger Green Safety Instruction Guide](#)
- [About Us](#)

RIFAK Vertical Ranger Green Safety Instruction Guide

Introduction

Thank you for choosing the RIFAK (RipAway IFAK), an individual first aid kit designed to provide essential medical supplies for the implementation of the MARCH protocol. This guide provides important safety information, usage instructions, and disposal guidelines to ensure safe and effective use of your RIFAK kit.

General Safety Guidelines

- **Product Safety:** The RIFAK is designed to ensure safe use and contains no hazardous materials. Always follow the instructions for safe handling and usage.
- **Enhanced Recalls:** In the event of a product recall, you will be notified through standardized notices. Please stay informed about any updates regarding your product.
- **Online Shopping Safety:** If purchased online, ensure that the retailer complies with all safety requirements.
- **Special Consumer Focus:** This product is not intended for use by children. Keep it out of reach and ensure that it is used only by trained adults.
- **EU Contact Point:** For safety inquiries, please refer to the manufacturer's contact information provided with the product.
- **Rapid Alerts:** Stay updated on any safety alerts related to your product through the EU Safety Gate platform.

Specific Safety Precautions for Use

- **Handling:** Always handle the RIFAK with care. Ensure that the kit is closed securely to prevent accidental opening.
- **Training:** Only individuals trained in first aid should use the medical supplies contained within the RIFAK.
- **Inspection:** Regularly inspect the kit for any damage or expired supplies. Replace any items as necessary.
- **Emergency Situations:** In case of an emergency, ensure that the RIFAK is easily accessible and that all users know how to operate the RipAway system.
- **Environmental Conditions:** Store the RIFAK in a cool, dry place, away from direct sunlight and extreme temperatures to maintain the integrity of the medical supplies.

Instructions for Installation and Usage

1. Installation:

- The RIFAK can be secured using the PALS/MOLLE system. Attach the kit to your gear using compatible straps.
- Ensure the kit is firmly attached and does not obstruct movement.

2. Usage:

- To access medical supplies, pull the RipAway handle. The internal insert will detach from the outer shell.
- In less than 6 seconds, you should be able to access all necessary medical supplies.
- Ensure that you are familiar with the layout of the kit and the location of each medical supply.

3. Additional Features:

- Utilize the external elastic cord to compress the size of the pocket or to carry a Tourniquet.
- You can attach additional pockets or tools (e.g., scissors, strapcutters) to the side PALS/MOLLE for enhanced functionality.

Disposal Instructions

- Dispose of the RIFAK according to local regulations for waste disposal.
- If the kit is damaged or expired, ensure that it is disposed of safely to prevent accidental use of the medical supplies.
- Recycle any components that are recyclable, in accordance with local recycling guidelines.

Contact Information for Further Support

For any further assistance or inquiries regarding the RIFAK Vertical Ranger Green, please refer to the contact information provided with your product packaging. It is important to reach out to the manufacturer for any safety concerns or questions regarding the use of the kit.

Thank you for your attention to these safety guidelines. Your safety and the effective use of the RIFAK kit are our top priorities.

About Us

Brownells Italia

Brownells Italia - The shop of the largest dealer of components and accessories for target shooting guns, defense and hunting

**Sede Legale e Operativa **

Brownells Italia Srl

Via Dante Alighieri 163

25068 Sarezzo Bs

P.Iva e C.F. IT02801420981

www.brownells.it